

Material Compatibility Status

Fresenius		
Equipment	Metrex Product	Compatibility Status
CRM III Outer Casing	CaviWipes1™	No information
CRM III Outer Casing	CaviWipes™ 2.0	No information
CRM III Outer Casing	CaviWipes™ HP	No information
CRM III Outer Casing	CaviWipes™ Bleach	No information
GranuFlo Concentrate Dissolution Unit II	CaviCide™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™	No information
GranuFlo Concentrate Dissolution Unit II	Envirocide™	No information
GranuFlo Concentrate Dissolution Unit II	CaviCide1™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes1™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ 2.0	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ HP	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ Bleach	Compatible
Liberty Cyclor	CaviCide™	Not recommended for use
Liberty Cyclor	CaviWipes™	Not recommended for use
Liberty Cyclor	Envirocide™	Not recommended for use
Liberty Cyclor	CaviCide1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Liberty Cyclor	CaviWipes1™	Not recommended for use
Liberty Cyclor	CaviWipes™ 2.0	Not recommended for use
Liberty Cyclor	CaviWipes™ HP	No information
Liberty Cyclor	CaviWipes™ Bleach	Not recommended for use
Newton Cyclor	CaviCide™	Not recommended for use
Newton Cyclor	CaviWipes™	Not recommended for use
Newton Cyclor	Envirocide™	Not recommended for use
Newton Cyclor	CaviCide1™	Not recommended for use
Newton Cyclor	CaviWipes1™	Not recommended for use
Newton Cyclor	CaviWipes™ 2.0	Not recommended for use
Newton Cyclor	CaviWipes™ HP	No information
Newton Cyclor	CaviWipes™ Bleach	Not recommended for use
2008K	CaviCide™	Not recommended for use
2008K	CaviWipes™	Not recommended for use
2008K	Envirocide™	Not recommended for use
2008K	CaviCide1™	Not recommended for use
2008K	CaviWipes1™	Not recommended for use
2008K	CaviWipes™ 2.0	Not recommended for use
2008K	CaviWipes™ HP	No information
2008K	CaviWipes™ Bleach	Not recommended for use
Crit-Line Monitor III Display	CaviCide™	No information
Crit-Line Monitor III Display	CaviWipes™	No information
Crit-Line Monitor III Display	Envirocide™	No information
Crit-Line Monitor III Display	CaviCide1™	No information
Crit-Line Monitor III Display	CaviWipes1™	No information
Crit-Line Monitor III Display	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Crit-Line Monitor III Display	CaviWipes™ HP	No information
Crit-Line Monitor III Display	CaviWipes™ Bleach	Not recommended for use
Crit-Line Monitor III Sensor clip	CaviCide™	Compatible
Crit-Line Monitor III Sensor clip	CaviWipes™	Compatible
Crit-Line Monitor III Sensor clip	Envirocide™	Compatible
Crit-Line Monitor III Sensor clip	CaviCide1™	No information
Crit-Line Monitor III Sensor clip	CaviWipes1™	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ 2.0	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ HP	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ Bleach	Compatible
CRM III Outer Casing	CaviCide™	Compatible
CRM III Outer Casing	CaviWipes™	Compatible
CRM III Outer Casing	Envirocide™	Compatible
CRM III Outer Casing	CaviCide1™	No information

Spacelabs		
Equipment	Metrex Product	Compatibility Status
Tri-Cuff	CaviCide™	No information
Tri-Cuff	CaviWipes™	No information
Tri-Cuff	Envirocide™	No information
Tri-Cuff	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Tri-Cuff	CaviWipes1™	No information
Tri-Cuff	CaviWipes™ 2.0	No information
Tri-Cuff	CaviWipes™ HP	No information
Tri-Cuff	CaviWipes™ Bleach	No information
TRU CUFF	CaviCide™	No information
TRU CUFF	CaviWipes™	No information
TRU CUFF	Envirocide™	No information
TRU CUFF	CaviCide1™	No information
TRU CUFF	CaviWipes1™	No information
TRU CUFF	CaviWipes™ 2.0	No information
TRU CUFF	CaviWipes™ HP	No information
TRU CUFF	CaviWipes™ Bleach	Compatible
TRU LINK	CaviCide™	No information
TRU LINK	CaviWipes™	No information
TRU LINK	Envirocide™	No information
TRU LINK	CaviCide1™	No information
TRU LINK	CaviWipes1™	No information
TRU LINK	CaviWipes™ 2.0	No information
TRU LINK	CaviWipes™ HP	No information
TRU LINK	CaviWipes™ Bleach	Compatible
Ultaview SL	CaviCide™	No information
Ultaview SL	CaviWipes™	No information
Ultaview SL	Envirocide™	No information
Ultaview SL	CaviCide1™	No information
Ultaview SL	CaviWipes1™	No information
Ultaview SL	CaviWipes™ 2.0	No information
Ultaview SL	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Ultaview SL	CaviWipes™ Bleach	Compatible
Ultraview Care Network	CaviCide™	No information
Ultraview Care Network	CaviWipes™	No information
Ultraview Care Network	Envirocide™	No information
Ultraview Care Network	CaviCide1™	No information
Ultraview Care Network	CaviWipes1™	No information
Ultraview Care Network	CaviWipes™ 2.0	No information
Ultraview Care Network	CaviWipes™ HP	No information
Ultraview Care Network	CaviWipes™ Bleach	No information
Ultraview Care Network/Monitor	CaviCide™	No information
Ultraview Care Network/Monitor	CaviWipes™	No information
Ultraview Care Network/Monitor	Envirocide™	No information
Ultraview Care Network/Monitor	CaviCide1™	No information
Ultraview Care Network/Monitor	CaviWipes1™	No information
Ultraview Care Network/Monitor	CaviWipes™ 2.0	No information
Ultraview Care Network/Monitor	CaviWipes™ HP	No information
Ultraview Care Network/Monitor	CaviWipes™ Bleach	No information
Ultraview Care Network/Transducers	CaviCide™	No information
Ultraview Care Network/Transducers	CaviWipes™	No information
Ultraview Care Network/Transducers	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Ultraview Care Network/Transducers	CaviCide1™	No information
Ultraview Care Network/Transducers	CaviWipes1™	No information
Ultraview Care Network/Transducers	CaviWipes™ 2.0	No information
Ultraview Care Network/Transducers	CaviWipes™ HP	No information
Ultraview Care Network/Transducers	CaviWipes™ Bleach	No information
Ultraview Care Network/Sensors	CaviCide™	No information
Ultraview Care Network/Sensors	CaviWipes™	No information
Ultraview Care Network/Sensors	Envirocide™	No information
Ultraview Care Network/Sensors	CaviCide1™	No information
Ultraview Care Network/Sensors	CaviWipes1™	No information
Ultraview Care Network/Sensors	CaviWipes™ 2.0	No information
Ultraview Care Network/Sensors	CaviWipes™ HP	No information
Ultraview Care Network/Sensors	CaviWipes™ Bleach	Compatible