

Material Compatibility Status

Dexis		
Equipment	Metrex Product	Compatibility Status
Gendex GXS-700	CaviCide™	Compatible
Gendex GXS-700	CaviWipes™	Compatible
Gendex GXS-700	Envirocide™	Compatible
Gendex GXS-700	CaviCide1™	No information
Gendex GXS-700	CaviWipes1™	No information
Gendex GXS-700	CaviWipes™ 2.0	No information
Gendex GXS-700	CaviWipes™ HP	No information
Gendex GXS-700	CaviWipes™ Bleach	No information

Geirger		
Equipment	Metrex Product	Compatibility Status
TCU 150	CaviCide™	No information
TCU 150	CaviWipes™	No information
TCU 150	Envirocide™	No information
TCU 150	CaviCide1™	No information
TCU 150	CaviWipes1™	No information
TCU 150	CaviWipes™ 2.0	No information
TCU 150	CaviWipes™ HP	No information
TCU 150	CaviWipes™ Bleach	No information

Posey		
-------	--	--

Equipment	Metrex Product	Compatibility Status
Bed 8060	CaviCide™	No information
Bed 8060	CaviWipes™	No information
Bed 8060	Envirocide™	No information
Bed 8060	CaviCide1™	No information
Bed 8060	CaviWipes1™	No information
Bed 8060	CaviWipes™ 2.0	No information
Bed 8060	CaviWipes™ HP	No information
Bed 8060	CaviWipes™ Bleach	No information
Bed 8070	CaviCide™	Compatible
Bed 8070	CaviWipes™	Compatible
Bed 8070	Envirocide™	Compatible
Bed 8070	CaviCide1™	No information
Bed 8070	CaviWipes1™	No information
Bed 8070	CaviWipes™ 2.0	Compatible
Bed 8070	CaviWipes™ HP	Compatible
Bed 8070	CaviWipes™ Bleach	Compatible
Floor cushions	CaviCide™	Compatible
Floor cushions	CaviWipes™	Compatible
Floor cushions	Envirocide™	Compatible
Floor cushions	CaviCide1™	Compatible
Floor cushions	CaviWipes1™	Compatible
Floor cushions	CaviWipes™ 2.0	Compatible
Floor cushions	CaviWipes™ HP	Compatible
Floor cushions	CaviWipes™ Bleach	Compatible