

Material Compatibility Status

Covidien/Medtronic/Puritan Bennett		
Equipment	Metrex Product	Compatibility Status
800 Series/840	CaviCide™	Compatible
800 Series/840	CaviWipes™	Compatible
800 Series/840	Envirocide™	Compatible
800 Series/840	CaviCide1™	No information
800 Series/840	CaviWipes1™	No information
800 Series/840	CaviWipes™ 2.0	No information
800 Series/840	CaviWipes™ HP	Compatible
800 Series/840	CaviWipes™ Bleach	Compatible

Wavesense		
Equipment	Metrex Product	Compatibility Status
Keynote Pro	CaviCide™	No information
Keynote Pro	CaviWipes™	No information
Keynote Pro	Envirocide™	No information
Keynote Pro	CaviCide1™	No information
Keynote Pro	CaviWipes1™	No information
Keynote Pro	CaviWipes™ 2.0	No information
Keynote Pro	CaviWipes™ HP	No information
Keynote Pro	CaviWipes™ Bleach	Compatible