

Material Compatibility Status

GE Healthcare		
Equipment	Metrex Product	Compatibility Status
RAB2-6-RS	CaviCide™	Compatible
RAB2-6-RS	CaviWipes™	Compatible
RAB2-6-RS	Envirocide™	Compatible
RAB2-6-RS	CaviCide1™	Compatible
RAB2-6-RS	CaviWipes1™	Compatible
RAB2-6-RS	CaviWipes™ 2.0	No information
RAB2-6-RS	CaviWipes™ HP	No information
RAB2-6-RS	CaviWipes™ Bleach	No information