

Material Compatibility Status

Abbott Healthcare		
Equipment	Metrex Product	Compatibility Status
Precision Xceed Pro	CaviCide™	Compatible
Precision Xceed Pro	CaviWipes™	Compatible
Precision Xceed Pro	Envirocide™	Compatible
Precision Xceed Pro	CaviCide1™	Not recommended for use
Precision Xceed Pro	CaviWipes1™	Not recommended for use
Precision Xceed Pro	CaviWipes™ 2.0	Compatible
Precision Xceed Pro	CaviWipes™ HP	Not recommended for use
Precision Xceed Pro	CaviWipes™ Bleach	Not recommended for use
Precision Xtra	CaviCide™	Not recommended for use
Precision Xtra	CaviWipes™	Not recommended for use
Precision Xtra	Envirocide™	Not recommended for use
Precision Xtra	CaviCide1™	Not recommended for use
Precision Xtra	CaviWipes1™	Not recommended for use
Precision Xtra	CaviWipes™ 2.0	Not recommended for use
Precision Xtra	CaviWipes™ HP	Not recommended for use
Precision Xtra	CaviWipes™ Bleach	Compatible
Flexible Quantum	CaviCide™	Not recommended for use
Flexible Quantum	CaviWipes™	Not recommended for use
Flexible Quantum	Envirocide™	Not recommended for use
Flexible Quantum	CaviCide1™	Not recommended for use
Flexible Quantum	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Flexible Quantum	CaviWipes™ 2.0	Not recommended for use
Flexible Quantum	CaviWipes™ HP	Not recommended for use
Flexible Quantum	CaviWipes™ Bleach	Compatible
Freestyle Flash	CaviCide™	Not recommended for use
Freestyle Flash	CaviWipes™	Not recommended for use
Freestyle Flash	Envirocide™	Not recommended for use
Freestyle Flash	CaviCide1™	Not recommended for use
Freestyle Flash	CaviWipes1™	Not recommended for use
Freestyle Flash	CaviWipes™ 2.0	Not recommended for use
Freestyle Flash	CaviWipes™ HP	Not recommended for use
Freestyle Flash	CaviWipes™ Bleach	Compatible
Freestyle Freedom	CaviCide™	Not recommended for use
Freestyle Freedom	CaviWipes™	Not recommended for use
Freestyle Freedom	Envirocide™	Not recommended for use
Freestyle Freedom	CaviCide1™	Not recommended for use
Freestyle Freedom	CaviWipes1™	Not recommended for use
Freestyle Freedom	CaviWipes™ 2.0	Not recommended for use
Freestyle Freedom	CaviWipes™ HP	Not recommended for use
Freestyle Freedom	CaviWipes™ Bleach	Compatible
Freestyle Lite	CaviWipes1™	Not recommended for use
Freestyle Lite	CaviWipes™ 2.0	Not recommended for use
Freestyle Lite	CaviWipes™ HP	Not recommended for use
Freestyle Lite	CaviWipes™ Bleach	Compatible
Optimum Xceed	CaviCide™	Not recommended for use
Optimum Xceed	CaviWipes™	Not recommended for use
Optimum Xceed	Envirocide™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Optimum Xceed	CaviCide1™	Not recommended for use
Optimum Xceed	CaviWipes1™	Not recommended for use
Optimum Xceed	CaviWipes™ 2.0	Not recommended for use
Optimum Xceed	CaviWipes™ HP	Not recommended for use
Optimum Xceed	CaviWipes™ Bleach	Compatible
Patrol ePump	CaviCide™	Not recommended for use
Patrol ePump	CaviWipes™	Not recommended for use
Patrol ePump	Envirocide™	Not recommended for use
Patrol ePump	CaviCide1™	Not recommended for use
Patrol ePump	CaviWipes1™	Not recommended for use
Patrol ePump	CaviWipes™ 2.0	Not recommended for use
Patrol ePump	CaviWipes™ HP	Not recommended for use
Patrol ePump	CaviWipes™ Bleach	Compatible
FreeStyle Freedom Lite	CaviCide™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™	Not recommended for use
FreeStyle Freedom Lite	Envirocide™	Not recommended for use
FreeStyle Freedom Lite	CaviCide1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ 2.0	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ HP	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ Bleach	Compatible
Freestyle Insulinx	CaviCide™	Not recommended for use
Freestyle Insulinx	CaviWipes™	Not recommended for use
Freestyle Insulinx	Envirocide™	Not recommended for use
Freestyle Insulinx	CaviCide1™	Not recommended for use
Freestyle Insulinx	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Freestyle Insulinx	CaviWipes™ 2.0	Not recommended for use
Freestyle Insulinx	CaviWipes™ HP	Not recommended for use
Freestyle Insulinx	CaviWipes™ Bleach	Not recommended for use
Freestyle Lite	CaviCide™	Not recommended for use
Freestyle Lite	CaviWipes™	Not recommended for use
Freestyle Lite	Envirocide™	Not recommended for use
Freestyle Lite	CaviCide1™	Not recommended for use

Cardiac science		
Equipment	Metrex Product	Compatibility Status
Power Heart AED G3 9300E	CaviCide™	No information
Power Heart AED G3 9300E	CaviWipes™	No information
Power Heart AED G3 9300E	Envirocide™	No information
Power Heart AED G3 9300E	CaviCide1™	No information
Power Heart AED G3 9300E	CaviWipes1™	No information
Power Heart AED G3 9300E	CaviWipes™ 2.0	No information
Power Heart AED G3 9300E	CaviWipes™ HP	No information
Power Heart AED G3 9300E	CaviWipes™ Bleach	No information
Quinton Eclipse Premier ECG	CaviCide™	No information
Quinton Eclipse Premier ECG	CaviWipes™	No information
Quinton Eclipse Premier ECG	Envirocide™	No information
Quinton Eclipse Premier ECG	CaviCide1™	No information
Quinton Eclipse Premier ECG	CaviWipes1™	No information
Quinton Eclipse Premier ECG	CaviWipes™ 2.0	No information
Quinton Eclipse Premier ECG	CaviWipes™ HP	No information
Quinton Eclipse Premier ECG	CaviWipes™ Bleach	No information