

# Material Compatibility Status

| Arizona Healthcare |                   |                      |
|--------------------|-------------------|----------------------|
| Equipment          | Metrex Product    | Compatibility Status |
| SpotOn Temp        | CaviCide™         | No information       |
| SpotOn Temp        | CaviWipes™        | No information       |
| SpotOn Temp        | Envirocide™       | No information       |
| SpotOn Temp        | CaviCide1™        | No information       |
| SpotOn Temp        | CaviWipes1™       | No information       |
| SpotOn Temp        | CaviWipes™ 2.0    | No information       |
| SpotOn Temp        | CaviWipes™ HP     | No information       |
| SpotOn Temp        | CaviWipes™ Bleach | No information       |

| Physio-Control |                   |                      |
|----------------|-------------------|----------------------|
| Equipment      | Metrex Product    | Compatibility Status |
| LifePak 12     | CaviCide™         | Compatible           |
| LifePak 12     | CaviWipes™        | Compatible           |
| LifePak 12     | Envirocide™       | Compatible           |
| LifePak 12     | CaviCide1™        | Compatible           |
| LifePak 12     | CaviWipes1™       | Compatible           |
| LifePak 12     | CaviWipes™ 2.0    | Compatible           |
| LifePak 12     | CaviWipes™ HP     | No information       |
| LifePak 12     | CaviWipes™ Bleach | No information       |
| LifePak 15     | CaviCide™         | Compatible           |
| LifePak 15     | CaviWipes™        | Compatible           |

| <b>Equipment</b>                          | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------------------------|-----------------------|-----------------------------|
| LifePak 15                                | Envirocide™           | Compatible                  |
| LifePak 15                                | CaviCide1™            | Compatible                  |
| LifePak 15                                | CaviWipes1™           | Compatible                  |
| LifePak 15                                | CaviWipes™ 2.0        | Compatible                  |
| LifePak 15                                | CaviWipes™ HP         | No information              |
| LifePak 15                                | CaviWipes™ Bleach     | No information              |
| LifePak CR Plus Case Only Disposable Pads | CaviCide™             | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes™            | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | Envirocide™           | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviCide1™            | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes1™           | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes™ 2.0        | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes™ HP         | No information              |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes™ Bleach     | No information              |
| LifePak Express Case only disposable pads | CaviCide™             | Compatible                  |
| LifePak Express Case only disposable pads | CaviWipes™            | Compatible                  |
| LifePak Express Case only disposable pads | Envirocide™           | Compatible                  |
| LifePak Express Case only disposable pads | CaviCide1™            | Compatible                  |

| <b>Equipment</b>                          | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------------------------|-----------------------|-----------------------------|
| LifePak Express Case only disposable pads | CaviWipes1™           | Compatible                  |
| LifePak Express Case only disposable pads | CaviWipes™ 2.0        | Compatible                  |
| LifePak Express Case only disposable pads | CaviWipes™ HP         | No information              |
| LifePak Express Case only disposable pads | CaviWipes™ Bleach     | No information              |