

Material Compatibility Status

Bear Medical Systems		
Equipment	Metrex Product	Compatibility Status
Bear Cub 750PSV/exterior	CaviCide™	No information
Bear Cub 750PSV/exterior	CaviWipes™	No information
Bear Cub 750PSV/exterior	Envirocide™	No information
Bear Cub 750PSV/exterior	CaviCide1™	No information
Bear Cub 750PSV/exterior	CaviWipes1™	No information
Bear Cub 750PSV/exterior	CaviWipes™ 2.0	No information
Bear Cub 750PSV/exterior	CaviWipes™ HP	No information
Bear Cub 750PSV/exterior	CaviWipes™ Bleach	No information
Bear Cub 750PSV/sensor	CaviCide™	No information
Bear Cub 750PSV/sensor	CaviWipes™	No information
Bear Cub 750PSV/sensor	Envirocide™	No information
Bear Cub 750PSV/sensor	CaviCide1™	No information
Bear Cub 750PSV/sensor	CaviWipes1™	No information
Bear Cub 750PSV/sensor	CaviWipes™ 2.0	No information
Bear Cub 750PSV/sensor	CaviWipes™ HP	No information
Bear Cub 750PSV/sensor	CaviWipes™ Bleach	No information

Carestream		
Equipment	Metrex Product	Compatibility Status
CS 3600 Family	CaviCide™	Compatible
CS 3600 Family	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
CS 3600 Family	Envirocide™	Compatible
CS 3600 Family	CaviCide1™	No information
CS 3700	CaviWipes1™	Compatible
CS 3700	CaviWipes™ 2.0	Compatible
CS 3700	CaviWipes™ HP	Compatible
CS 3700	CaviWipes™ Bleach	Compatible
CS 8100	CaviCide™	Compatible
CS 8100	CaviWipes™	Compatible
CS 8100	Envirocide™	Compatible
CS 8100	CaviCide1™	Compatible
CS 8100	CaviWipes1™	Compatible
CS 8100	CaviWipes™ 2.0	Compatible
CS 8100	CaviWipes™ HP	Compatible
CS 8100	CaviWipes™ Bleach	Compatible
CS 9100	CaviCide™	No information
CS 9100	CaviWipes™	No information
CS 9100	Envirocide™	No information
CS 9100	CaviCide1™	No information
CS 9100	CaviWipes1™	No information
CS 9100	CaviWipes™ 2.0	No information
CS 9100	CaviWipes™ HP	No information
CS 9100	CaviWipes™ Bleach	No information
CS 9600	CaviCide™	Compatible
CS 9600	CaviWipes™	Compatible
CS 9600	Envirocide™	Compatible
CS 9600	CaviCide1™	Compatible

Equipment	Metrex Product	Compatibility Status
CS 9600	CaviWipes1™	Compatible
CS 9600	CaviWipes™ 2.0	Compatible
CS 9600	CaviWipes™ HP	Compatible
CS 9600	CaviWipes™ Bleach	Compatible
RVG Series IntraOral Scanner	CaviCide™	Compatible
RVG Series IntraOral Scanner	CaviWipes™	Compatible
RVG Series IntraOral Scanner	Envirocide™	Compatible
RVG Series IntraOral Scanner	CaviCide1™	Compatible
RVG Series IntraOral Scanner	CaviWipes1™	Compatible
RVG Series IntraOral Scanner	CaviWipes™ 2.0	Compatible
RVG Series IntraOral Scanner	CaviWipes™ HP	Compatible
RVG Series IntraOral Scanner	CaviWipes™ Bleach	Compatible
DRX Equipment	CaviCide™	No information
DRX Equipment	CaviWipes™	No information
DRX Equipment	Envirocide™	No information
DRX Equipment	CaviCide1™	Compatible
DRX Equipment	CaviWipes1™	Compatible
DRX Equipment	CaviWipes™ 2.0	No information
DRX Equipment	CaviWipes™ HP	Compatible
DRX Equipment	CaviWipes™ Bleach	Compatible
DryView Exterior	CaviCide™	No information
DryView Exterior	CaviWipes™	No information
DryView Exterior	Envirocide™	No information
DryView Exterior	CaviCide1™	Compatible
DryView Exterior	CaviWipes1™	Compatible
DryView Exterior	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
DryView Exterior	CaviWipes™ HP	Compatible
DryView Exterior	CaviWipes™ Bleach	Compatible
DRX Revolution	CaviCide™	No information
DRX Revolution	CaviWipes™	No information
DRX Revolution	Envirocide™	No information
DRX Revolution	CaviCide1™	Compatible
DRX Revolution	CaviWipes1™	Compatible
DRX Revolution	CaviWipes™ 2.0	No information
DRX Revolution	CaviWipes™ HP	No information
DRX Revolution	CaviWipes™ Bleach	No information
DRX-1	CaviCide™	No information
DRX-1	CaviWipes™	No information
DRX-1	Envirocide™	No information
DRX-1	CaviCide1™	Compatible
DRX-1	CaviWipes1™	Compatible
DRX-1	CaviWipes™ 2.0	No information
DRX-1	CaviWipes™ HP	Compatible
DRX-1	CaviWipes™ Bleach	Compatible
CS 1200	CaviCide™	No information
CS 1200	CaviWipes™	No information
CS 1200	Envirocide™	No information
CS 1200	CaviCide1™	No information
CS 1200	CaviWipes1™	No information
CS 1200	CaviWipes™ 2.0	No information
CS 1200	CaviWipes™ HP	No information
CS 1200	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
CS 3600 Family	CaviWipes1™	No information
CS 3600 Family	CaviWipes™ 2.0	No information
CS 3600 Family	CaviWipes™ HP	No information
CS 3600 Family	CaviWipes™ Bleach	Compatible
CS 3800	CaviCide™	Compatible
CS 3800	CaviWipes™	Compatible
CS 3800	Envirocide™	Compatible
CS 3800	CaviCide1™	No information
CS 3800	CaviWipes1™	No information
CS 3800	CaviWipes™ 2.0	No information
CS 3800	CaviWipes™ HP	No information
CS 3800	CaviWipes™ Bleach	No information
DR	CaviCide™	No information
DR	CaviWipes™	No information
DR	Envirocide™	No information
DR	CaviCide1™	Compatible
DR	CaviWipes1™	Compatible
DR	CaviWipes™ 2.0	No information
DR	CaviWipes™ HP	Compatible
DR	CaviWipes™ Bleach	Compatible
CS 1500	CaviCide™	No information
CS 1500	CaviWipes™	No information
CS 1500	Envirocide™	No information
CS 1500	CaviCide1™	No information
CS 1500	CaviWipes1™	No information
CS 1500	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
CS 1500	CaviWipes™ HP	No information
CS 1500	CaviWipes™ Bleach	Compatible
CS 3600	CaviCide™	Compatible
CS 3600	CaviWipes™	Compatible
CS 3600	Envirocide™	Compatible
CS 3600	CaviCide1™	Compatible
CS 3600	CaviWipes1™	Compatible
CS 3600	CaviWipes™ 2.0	Compatible
CS 3600	CaviWipes™ HP	Compatible
CS 3600	CaviWipes™ Bleach	Compatible
CS 3700	CaviCide™	Compatible
CS 3700	CaviWipes™	Compatible
CS 3700	Envirocide™	Compatible
CS 3700	CaviCide1™	Compatible

Hospira		
Equipment	Metrex Product	Compatibility Status
Plum A+3	CaviCide™	No information
Plum A+3	CaviWipes™	No information
Plum A+3	Envirocide™	No information
Plum A+3	CaviCide1™	No information
Plum A+3	CaviWipes1™	No information
Plum A+3	CaviWipes™ 2.0	No information
Plum A+3	CaviWipes™ HP	No information
Plum A+3	CaviWipes™ Bleach	Compatible
Plum XL3	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Plum XL3	CaviWipes™	No information
Plum XL3	Envirocide™	No information
Plum XL3	CaviCide1™	No information
Plum XL3	CaviWipes1™	No information
Plum XL3	CaviWipes™ 2.0	No information
Plum XL3	CaviWipes™ HP	No information
Plum XL3	CaviWipes™ Bleach	Compatible
Plum XL3 Micro/Macro	CaviCide™	No information
Plum XL3 Micro/Macro	CaviWipes™	No information
Plum XL3 Micro/Macro	Envirocide™	No information
Plum XL3 Micro/Macro	CaviCide1™	No information
Plum XL3 Micro/Macro	CaviWipes1™	No information
Plum XL3 Micro/Macro	CaviWipes™ 2.0	No information
Plum XL3 Micro/Macro	CaviWipes™ HP	No information
Plum XL3 Micro/Macro	CaviWipes™ Bleach	Compatible
Sapphire Pump	CaviCide™	No information
Sapphire Pump	CaviWipes™	No information
Sapphire Pump	Envirocide™	No information
Sapphire Pump	CaviCide1™	No information
Sapphire Pump	CaviWipes1™	No information
Sapphire Pump	CaviWipes™ 2.0	No information
Sapphire Pump	CaviWipes™ HP	No information
Sapphire Pump	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027-04-79 & above	CaviCide™	No information
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027-04-79 & above	CaviWipes™	No information
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027-04-79 & above	Envirocide™	No information
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027-04-79 & above	CaviCide1™	No information
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027-04-79 & above	CaviWipes1™	No information
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027-04-79 & above	CaviWipes™ 2.0	No information
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027-04-79 & above	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Symbiq One Channel Infuser List number 16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79 & above	CaviWipes™ Bleach	Compatible
LifeCare PCA Llist 2709-04	CaviCide™	No information
LifeCare PCA Llist 2709-04	CaviWipes™	No information
LifeCare PCA Llist 2709-04	Envirocide™	No information
LifeCare PCA Llist 2709-04	CaviCide1™	No information
LifeCare PCA Llist 2709-04	CaviWipes1™	No information
LifeCare PCA Llist 2709-04	CaviWipes™ 2.0	No information
LifeCare PCA Llist 2709-04	CaviWipes™ HP	No information
LifeCare PCA Llist 2709-04	CaviWipes™ Bleach	Compatible
Omni Flow 4000Plus	CaviCide™	No information
Omni Flow 4000Plus	CaviWipes™	No information
Omni Flow 4000Plus	Envirocide™	No information
Omni Flow 4000Plus	CaviCide1™	No information
Omni Flow 4000Plus	CaviWipes1™	No information
Omni Flow 4000Plus	CaviWipes™ 2.0	No information
Omni Flow 4000Plus	CaviWipes™ HP	No information
Omni Flow 4000Plus	CaviWipes™ Bleach	Compatible
Plum 360	CaviCide™	No information
Plum 360	CaviWipes™	No information
Plum 360	Envirocide™	No information
Plum 360	CaviCide1™	No information
Plum 360	CaviWipes1™	No information
Plum 360	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Plum 360	CaviWipes™ HP	No information
Plum 360	CaviWipes™ Bleach	Compatible
Plum A+	CaviCide™	No information
Plum A+	CaviWipes™	No information
Plum A+	Envirocide™	No information
Plum A+	CaviCide1™	No information
Plum A+	CaviWipes1™	No information
Plum A+	CaviWipes™ 2.0	No information
Plum A+	CaviWipes™ HP	No information
Plum A+	CaviWipes™ Bleach	Compatible
Plum A+ W/Hospira MedNet	CaviCide™	No information
Plum A+ W/Hospira MedNet	CaviWipes™	No information
Plum A+ W/Hospira MedNet	Envirocide™	No information
Plum A+ W/Hospira MedNet	CaviCide1™	No information
Plum A+ W/Hospira MedNet	CaviWipes1™	No information
Plum A+ W/Hospira MedNet	CaviWipes™ 2.0	No information
Plum A+ W/Hospira MedNet	CaviWipes™ HP	No information
Plum A+ W/Hospira MedNet	CaviWipes™ Bleach	Compatible
Acclaim Encore	CaviCide™	No information
Acclaim Encore	CaviWipes™	No information
Acclaim Encore	Envirocide™	No information
Acclaim Encore	CaviCide1™	No information
Acclaim Encore	CaviWipes1™	No information
Acclaim Encore	CaviWipes™ 2.0	No information
Acclaim Encore	CaviWipes™ HP	No information
Acclaim Encore	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
GemStar	CaviCide™	No information
GemStar	CaviWipes™	No information
GemStar	Envirocide™	No information
GemStar	CaviCide1™	No information
GemStar	CaviWipes1™	No information
GemStar	CaviWipes™ 2.0	No information
GemStar	CaviWipes™ HP	No information
GemStar	CaviWipes™ Bleach	Compatible

Nonin		
Equipment	Metrex Product	Compatibility Status
9843 Oximeter	CaviCide™	No information
9843 Oximeter	CaviWipes™	No information
9843 Oximeter	Envirocide™	No information
9843 Oximeter	CaviCide1™	No information
9843 Oximeter	CaviWipes1™	No information
9843 Oximeter	CaviWipes™ 2.0	No information
9843 Oximeter	CaviWipes™ HP	No information
9843 Oximeter	CaviWipes™ Bleach	No information