

Material Compatibility Status

Honeywell		
Equipment	Metrex Product	Compatibility Status
Voyager XP 1472g	CaviCide™	No information
Voyager XP 1472g	CaviWipes™	No information
Voyager XP 1472g	Envirocide™	No information
Voyager XP 1472g	CaviCide1™	No information
Voyager XP 1472g	CaviWipes1™	No information
Voyager XP 1472g	CaviWipes™ 2.0	No information
Voyager XP 1472g	CaviWipes™ HP	Compatible
Voyager XP 1472g	CaviWipes™ Bleach	No information
Xenon 1950g.bf	CaviCide™	No information
Xenon 1950g.bf	CaviWipes™	No information
Xenon 1950g.bf	Envirocide™	No information
Xenon 1950g.bf	CaviCide1™	No information
Xenon 1950g.bf	CaviWipes1™	No information
Xenon 1950g.bf	CaviWipes™ 2.0	No information
Xenon 1950g.bf	CaviWipes™ HP	Compatible
Xenon 1950g.bf	CaviWipes™ Bleach	No information
Xenon XP 1950g	CaviCide™	Compatible
Xenon XP 1950g	CaviWipes™	Compatible
Xenon XP 1950g	Envirocide™	Compatible
Xenon XP 1950g	CaviCide1™	No information
Xenon XP 1950g	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Xenon XP 1950g	CaviWipes™ 2.0	No information
Xenon XP 1950g	CaviWipes™ HP	No information
Xenon XP 1950g	CaviWipes™ Bleach	Compatible
Voyager XP 1470g	CaviCide™	No information
Voyager XP 1470g	CaviWipes™	No information
Voyager XP 1470g	Envirocide™	No information
Voyager XP 1470g	CaviCide1™	No information
Voyager XP 1470g	CaviWipes1™	No information
Voyager XP 1470g	CaviWipes™ 2.0	No information
Voyager XP 1470g	CaviWipes™ HP	Compatible
Voyager XP 1470g	CaviWipes™ Bleach	No information

WhiteHall MFG		
Equipment	Metrex Product	Compatibility Status
P-10-M-SDP	CaviCide™	No information
P-10-M-SDP	CaviWipes™	No information
P-10-M-SDP	Envirocide™	No information
P-10-M-SDP	CaviCide1™	No information
P-10-M-SDP	CaviWipes1™	No information
P-10-M-SDP	CaviWipes™ 2.0	No information
P-10-M-SDP	CaviWipes™ HP	No information
P-10-M-SDP	CaviWipes™ Bleach	No information