

Material Compatibility Status

MTI		
Equipment	Metrex Product	Compatibility Status
401	CaviWipes1™	No information
401	CaviWipes™ 2.0	No information
401	CaviWipes™ HP	No information
401	CaviWipes™ Bleach	Compatible
401 chairs	CaviCide™	No information
401 chairs	CaviWipes™	No information
401 chairs	Envirocide™	No information
401 chairs	CaviCide1™	No information
401 chairs	CaviWipes1™	No information
401 chairs	CaviWipes™ 2.0	No information
401 chairs	CaviWipes™ HP	No information
401 chairs	CaviWipes™ Bleach	Compatible
401	CaviCide™	No information
401	CaviWipes™	No information
401	Envirocide™	No information
401	CaviCide1™	No information