

Material Compatibility Status

NDS		
Equipment	Metrex Product	Compatibility Status
Radiance 19" Medical Display	CaviCide™	Compatible
Radiance 19" Medical Display	CaviWipes™	Compatible
Radiance 19" Medical Display	Envirocide™	Compatible
Radiance 19" Medical Display	CaviCide1™	Not recommended for use
Radiance 19" Medical Display	CaviWipes1™	Not recommended for use
Radiance 19" Medical Display	CaviWipes™ 2.0	Not recommended for use
Radiance 19" Medical Display	CaviWipes™ HP	Compatible
Radiance 19" Medical Display	CaviWipes™ Bleach	Compatible