

Material Compatibility Status

Ohio Medical		
Equipment	Metrex Product	Compatibility Status
CARE-E-VAC 3	CaviWipes1™	No information
CARE-E-VAC 3	CaviWipes™ 2.0	No information
CARE-E-VAC 3	CaviWipes™ HP	No information
CARE-E-VAC 3	CaviWipes™ Bleach	No information
Push to Set Intermittent Suction Unit	CaviCide™	Compatible
Push to Set Intermittent Suction Unit	CaviWipes™	Compatible
Push to Set Intermittent Suction Unit	Envirocide™	Compatible
Push to Set Intermittent Suction Unit	CaviCide1™	No information
Push to Set Intermittent Suction Unit	CaviWipes1™	No information
Push to Set Intermittent Suction Unit	CaviWipes™ 2.0	No information
Push to Set Intermittent Suction Unit	CaviWipes™ HP	Compatible
Push to Set Intermittent Suction Unit	CaviWipes™ Bleach	Compatible
CARE-E-VAC 3	CaviCide™	No information
CARE-E-VAC 3	CaviWipes™	No information
CARE-E-VAC 3	Envirocide™	No information
CARE-E-VAC 3	CaviCide1™	No information