

Material Compatibility Status

Posey		
Equipment	Metrex Product	Compatibility Status
Bed 8060	CaviCide™	No information
Bed 8070	CaviCide™	Compatible
Floor cushions	CaviCide™	Compatible

Sirona		
Equipment	Metrex Product	Compatibility Status
Schick Sensor	CaviCide™	Compatible
XIOS XG sensor & cable	CaviCide™	Compatible
Orthophos SL 2D	CaviCide™	Compatible
Orthophos SL 3D	CaviCide™	Compatible