

Material Compatibility Status

Align Technology		
Equipment	Metrex Product	Compatibility Status
iTero element	CaviCide™	Compatible
iTero element	CaviWipes™	Compatible
iTero element	Envirocide™	Compatible
iTero element	CaviCide1™	Compatible
iTero element	CaviWipes1™	Compatible
iTero element	CaviWipes™ 2.0	No information
iTero element	CaviWipes™ HP	No information
iTero element	CaviWipes™ Bleach	No information

B Medical systems		
Equipment	Metrex Product	Compatibility Status
All models	CaviCide™	No information
All models	CaviWipes™	No information
All models	Envirocide™	No information
All models	CaviCide1™	No information
All models	CaviWipes1™	No information
All models	CaviWipes™ 2.0	No information
All models	CaviWipes™ HP	No information
All models	CaviWipes™ Bleach	No information

WelchAllyn/Hill-Rom		
---------------------	--	--

Equipment	Metrex Product	Compatibility Status
GS 777 Wall Transformer	CaviCide™	Compatible
GS 777 Wall Transformer	CaviWipes™	Compatible
GS 777 Wall Transformer	Envirocide™	Compatible
GS 777 Wall Transformer	CaviCide1™	No information
GS 777 Wall Transformer	CaviWipes1™	No information
GS 777 Wall Transformer	CaviWipes™ 2.0	No information
GS 777 Wall Transformer	CaviWipes™ HP	No information
GS 777 Wall Transformer	CaviWipes™ Bleach	Compatible
Model 71110	CaviCide™	No information
Model 71110	CaviWipes™	No information
Model 71110	Envirocide™	No information
Model 71110	CaviCide1™	No information
Model 71110	CaviWipes1™	No information
Model 71110	CaviWipes™ 2.0	No information
Model 71110	CaviWipes™ HP	No information
Model 71110	CaviWipes™ Bleach	Compatible
Model 7114X	CaviCide™	No information
Model 7114X	CaviWipes™	No information
Model 7114X	Envirocide™	No information
Model 7114X	CaviCide1™	No information
Model 7114X	CaviWipes1™	No information
Model 7114X	CaviWipes™ 2.0	No information
Model 7114X	CaviWipes™ HP	No information
Model 7114X	CaviWipes™ Bleach	No information
Wireless Acquisition Module UTK	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Wireless Acquisition Module UTK	CaviWipes™	No information
Wireless Acquisition Module UTK	Envirocide™	No information
Wireless Acquisition Module UTK	CaviCide1™	No information
Wireless Acquisition Module UTK	CaviWipes1™	No information
Wireless Acquisition Module UTK	CaviWipes™ 2.0	No information
Wireless Acquisition Module UTK	CaviWipes™ HP	No information
Wireless Acquisition Module UTK	CaviWipes™ Bleach	Compatible