

Material Compatibility Status

AgaMatrix		
Equipment	Metrex Product	Compatibility Status
Jazz	CaviCide™	Not recommended for use
Jazz	CaviWipes™	Not recommended for use
Jazz	Envirocide™	Not recommended for use
Jazz	CaviCide1™	Not recommended for use
Jazz	CaviWipes1™	Not recommended for use
Jazz	CaviWipes™ 2.0	Not recommended for use
Jazz	CaviWipes™ HP	Not recommended for use
Jazz	CaviWipes™ Bleach	Not recommended for use

BD		
Equipment	Metrex Product	Compatibility Status
Pyxis Supply Technologies	CaviCide™	Compatible
Pyxis Supply Technologies	CaviWipes™	Compatible
Pyxis Supply Technologies	Envirocide™	Compatible
Pyxis Supply Technologies	CaviCide1™	No information
Pyxis Supply Technologies	CaviWipes1™	No information
Pyxis Supply Technologies	CaviWipes™ 2.0	No information
Pyxis Supply Technologies	CaviWipes™ HP	No information
Pyxis Supply Technologies	CaviWipes™ Bleach	No information
Alaris System	CaviCide™	No information
Alaris System	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Alaris System	Envirocide™	No information
Alaris System	CaviCide1™	No information
Alaris System	CaviWipes1™	No information
Alaris System	CaviWipes™ 2.0	No information
Alaris System	CaviWipes™ HP	No information
Alaris System	CaviWipes™ Bleach	No information
Pyxis Anesthesia System	CaviCide™	No information
Pyxis Anesthesia System	CaviWipes™	No information
Pyxis Anesthesia System	Envirocide™	No information
Pyxis Anesthesia System	CaviCide1™	No information
Pyxis Anesthesia System	CaviWipes1™	No information
Pyxis Anesthesia System	CaviWipes™ 2.0	No information
Pyxis Anesthesia System	CaviWipes™ HP	No information
Pyxis Anesthesia System	CaviWipes™ Bleach	No information
Pyxis MedStation	CaviCide™	No information
Pyxis MedStation	CaviWipes™	No information
Pyxis MedStation	Envirocide™	No information
Pyxis MedStation	CaviCide1™	No information
Pyxis MedStation	CaviWipes1™	No information
Pyxis MedStation	CaviWipes™ 2.0	No information
Pyxis MedStation	CaviWipes™ HP	No information
Pyxis MedStation	CaviWipes™ Bleach	No information
Pyxis Supply Technologies	CaviCide™	Compatible
Pyxis Supply Technologies	CaviWipes™	Compatible
Pyxis Supply Technologies	Envirocide™	Compatible
Pyxis Supply Technologies	CaviCide1™	No information
Pyxis Supply Technologies	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Pyxis Supply Technologies	CaviWipes™ 2.0	No information
Pyxis Supply Technologies	CaviWipes™ HP	No information
Pyxis Supply Technologies	CaviWipes™ Bleach	No information