

Material Compatibility Status

CIVCO Medical Solutions		
Equipment	Metrex Product	Compatibility Status
Astra VR	CaviWipes1™	No information
Astra VR	CaviWipes™ 2.0	No information
Astra VR	CaviWipes™ HP	No information
Astra VR	CaviWipes™ Bleach	No information
Astra VR Ultra	CaviCide™	No information
Astra VR Ultra	CaviWipes™	No information
Astra VR Ultra	Envirocide™	No information
Astra VR Ultra	CaviCide1™	No information
Astra VR Ultra	CaviWipes1™	No information
Astra VR Ultra	CaviWipes™ 2.0	No information
Astra VR Ultra	CaviWipes™ HP	No information
Astra VR Ultra	CaviWipes™ Bleach	No information
Astra TEE	CaviCide™	No information
Astra TEE	CaviWipes™	No information
Astra TEE	Envirocide™	No information
Astra TEE	CaviCide1™	No information
Astra TEE	CaviWipes1™	No information
Astra TEE	CaviWipes™ 2.0	No information
Astra TEE	CaviWipes™ HP	No information
Astra TEE	CaviWipes™ Bleach	No information
Astra TEE Ultra	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Astra TEE Ultra	CaviWipes™	No information
Astra TEE Ultra	Envirocide™	No information
Astra TEE Ultra	CaviCide1™	No information
Astra TEE Ultra	CaviWipes1™	No information
Astra TEE Ultra	CaviWipes™ 2.0	No information
Astra TEE Ultra	CaviWipes™ HP	No information
Astra TEE Ultra	CaviWipes™ Bleach	No information
Astra VR	CaviCide™	No information
Astra VR	CaviWipes™	No information
Astra VR	Envirocide™	No information
Astra VR	CaviCide1™	No information

Posey		
Equipment	Metrex Product	Compatibility Status
Bed 8060	CaviCide™	No information
Bed 8060	CaviWipes™	No information
Bed 8060	Envirocide™	No information
Bed 8060	CaviCide1™	No information
Bed 8060	CaviWipes1™	No information
Bed 8060	CaviWipes™ 2.0	No information
Bed 8060	CaviWipes™ HP	No information
Bed 8060	CaviWipes™ Bleach	No information
Bed 8070	CaviCide™	Compatible
Bed 8070	CaviWipes™	Compatible
Bed 8070	Envirocide™	Compatible
Bed 8070	CaviCide1™	No information
Bed 8070	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Bed 8070	CaviWipes™ 2.0	Compatible
Bed 8070	CaviWipes™ HP	Compatible
Bed 8070	CaviWipes™ Bleach	Compatible
Floor cushions	CaviCide™	Compatible
Floor cushions	CaviWipes™	Compatible
Floor cushions	Envirocide™	Compatible
Floor cushions	CaviCide1™	Compatible
Floor cushions	CaviWipes1™	Compatible
Floor cushions	CaviWipes™ 2.0	Compatible
Floor cushions	CaviWipes™ HP	Compatible
Floor cushions	CaviWipes™ Bleach	Compatible