

Material Compatibility Status

Intuitive Surgical		
Equipment	Metrex Product	Compatibility Status
da Vinci Aecessories - Energy activation cable	CaviCide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Aecessories - Energy activation cable	Envirocide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviCide1™	No information
da Vinci Aecessories - Energy activation cable	CaviWipes1™	No information
da Vinci Aecessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Aecessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Aecessories - Energy activation cable	CaviWipes™ Bleach	No information
da Vinci Aecessories - Energy activation cable	CaviCide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Aecessories - Energy activation cable	Envirocide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
daVinci S System	CaviCide™	Compatible
daVinci S System	CaviWipes™	Compatible
daVinci S System	Envirocide™	Compatible
daVinci S System	CaviCide1™	No information
daVinci S System	CaviWipes1™	No information
daVinci S System	CaviWipes™ 2.0	No information
daVinci S System	CaviWipes™ HP	No information
daVinci S System	CaviWipes™ Bleach	No information
daVinci Xi System	CaviCide™	Compatible
daVinci Xi System	CaviWipes™	Compatible
daVinci Xi System	Envirocide™	Compatible
daVinci Xi System	CaviCide1™	No information
daVinci Xi System	CaviWipes1™	No information
daVinci Xi System	CaviWipes™ 2.0	No information
daVinci Xi System	CaviWipes™ HP	No information
daVinci Xi System	CaviWipes™ Bleach	No information

WelchAllyn/Hill-Rom		
Equipment	Metrex Product	Compatibility Status
GS 777 Wall Transformer	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
GS 777 Wall Transformer	CaviWipes™	Compatible
GS 777 Wall Transformer	Envirocide™	Compatible
GS 777 Wall Transformer	CaviCide1™	No information
GS 777 Wall Transformer	CaviWipes1™	No information
GS 777 Wall Transformer	CaviWipes™ 2.0	No information
GS 777 Wall Transformer	CaviWipes™ HP	No information
GS 777 Wall Transformer	CaviWipes™ Bleach	Compatible
Model 71110	CaviCide™	No information
Model 71110	CaviWipes™	No information
Model 71110	Envirocide™	No information
Model 71110	CaviCide1™	No information
Model 71110	CaviWipes1™	No information
Model 71110	CaviWipes™ 2.0	No information
Model 71110	CaviWipes™ HP	No information
Model 71110	CaviWipes™ Bleach	Compatible
Model 7114X	CaviCide™	No information
Model 7114X	CaviWipes™	No information
Model 7114X	Envirocide™	No information
Model 7114X	CaviCide1™	No information
Model 7114X	CaviWipes1™	No information
Model 7114X	CaviWipes™ 2.0	No information
Model 7114X	CaviWipes™ HP	No information
Model 7114X	CaviWipes™ Bleach	No information
Wireless Acquisition Module UTK	CaviCide™	No information
Wireless Acquisition Module UTK	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Wireless Acquisition Module UTK	Envirocide™	No information
Wireless Acquisition Module UTK	CaviCide1™	No information
Wireless Acquisition Module UTK	CaviWipes1™	No information
Wireless Acquisition Module UTK	CaviWipes™ 2.0	No information
Wireless Acquisition Module UTK	CaviWipes™ HP	No information
Wireless Acquisition Module UTK	CaviWipes™ Bleach	Compatible