

Material Compatibility Status

Intuitive Surgical		
Equipment	Metrex Product	Compatibility Status
da Vinci Aecessories - Energy activation cable	CaviCide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Aecessories - Energy activation cable	Envirocide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviCide1™	No information
da Vinci Aecessories - Energy activation cable	CaviWipes1™	No information
da Vinci Aecessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Aecessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Aecessories - Energy activation cable	CaviWipes™ Bleach	No information
da Vinci Aecessories - Energy activation cable	CaviCide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Aecessories - Energy activation cable	Envirocide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
daVinci S System	CaviCide™	Compatible
daVinci S System	CaviWipes™	Compatible
daVinci S System	Envirocide™	Compatible
daVinci S System	CaviCide1™	No information
daVinci S System	CaviWipes1™	No information
daVinci S System	CaviWipes™ 2.0	No information
daVinci S System	CaviWipes™ HP	No information
daVinci S System	CaviWipes™ Bleach	No information
daVinci Xi System	CaviCide™	Compatible
daVinci Xi System	CaviWipes™	Compatible
daVinci Xi System	Envirocide™	Compatible
daVinci Xi System	CaviCide1™	No information
daVinci Xi System	CaviWipes1™	No information
daVinci Xi System	CaviWipes™ 2.0	No information
daVinci Xi System	CaviWipes™ HP	No information
daVinci Xi System	CaviWipes™ Bleach	No information

WhiteHall MFG		
Equipment	Metrex Product	Compatibility Status
P-10-M-SDP	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
P-10-M-SDP	CaviWipes™	No information
P-10-M-SDP	Envirocide™	No information
P-10-M-SDP	CaviCide1™	No information
P-10-M-SDP	CaviWipes1™	No information
P-10-M-SDP	CaviWipes™ 2.0	No information
P-10-M-SDP	CaviWipes™ HP	No information
P-10-M-SDP	CaviWipes™ Bleach	No information