

Material Compatibility Status

Covaris		
Equipment	Metrex Product	Compatibility Status
cryoPREP Extraction system	CaviCide™	No information
cryoPREP Extraction system	CaviWipes™	No information
cryoPREP Extraction system	Envirocide™	No information
cryoPREP Extraction system	CaviCide1™	No information
cryoPREP Extraction system	CaviWipes1™	No information
cryoPREP Extraction system	CaviWipes™ 2.0	No information
cryoPREP Extraction system	CaviWipes™ HP	No information
cryoPREP Extraction system	CaviWipes™ Bleach	Compatible

Intuitive Surgical		
Equipment	Metrex Product	Compatibility Status
da Vinci Aecessories - Energy activation cable	CaviCide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Aecessories - Energy activation cable	Envirocide™	Compatible
da Vinci Aecessories - Energy activation cable	CaviCide1™	No information
da Vinci Aecessories - Energy activation cable	CaviWipes1™	No information
da Vinci Aecessories - Energy activation cable	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
da Vinci Acessories - Energy activation cable	CaviCide™	Compatible
da Vinci Acessories - Energy activation cable	CaviWipes™	Compatible
da Vinci Acessories - Energy activation cable	Envirocide™	Compatible
da Vinci Acessories - Energy activation cable	CaviCide1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes1™	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ 2.0	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ HP	No information
da Vinci Acessories - Energy activation cable	CaviWipes™ Bleach	No information
daVinci S System	CaviCide™	Compatible
daVinci S System	CaviWipes™	Compatible
daVinci S System	Envirocide™	Compatible
daVinci S System	CaviCide1™	No information
daVinci S System	CaviWipes1™	No information
daVinci S System	CaviWipes™ 2.0	No information
daVinci S System	CaviWipes™ HP	No information
daVinci S System	CaviWipes™ Bleach	No information
daVinci Xi System	CaviCide™	Compatible
daVinci Xi System	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
daVinci Xi System	Envirocide™	Compatible
daVinci Xi System	CaviCide1™	No information
daVinci Xi System	CaviWipes1™	No information
daVinci Xi System	CaviWipes™ 2.0	No information
daVinci Xi System	CaviWipes™ HP	No information
daVinci Xi System	CaviWipes™ Bleach	No information