

Material Compatibility Status

Accuvein		
Equipment	Metrex Product	Compatibility Status
AV400	CaviWipes1™	Not recommended for use
AV400	CaviWipes™ 2.0	No information
AV400	CaviWipes™ HP	No information
AV400	CaviWipes™ Bleach	Compatible
AV400	CaviCide™	No information
AV400	CaviWipes™	No information
AV400	Envirocide™	No information
AV400	CaviCide1™	No information
AV400	CaviWipes1™	No information
AV400	CaviWipes™ 2.0	No information
AV400	CaviWipes™ HP	No information
AV400	CaviWipes™ Bleach	No information
AV 300	CaviCide™	No information
AV 300	CaviWipes™	No information
AV 300	Envirocide™	No information
AV 300	CaviCide1™	Not recommended for use
AV 300	CaviWipes1™	Not recommended for use
AV 300	CaviWipes™ 2.0	No information
AV 300	CaviWipes™ HP	No information
AV 300	CaviWipes™ Bleach	Compatible
AV300	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
AV300	CaviWipes™	No information
AV300	Envirocide™	No information
AV300	CaviCide1™	No information
AV300	CaviWipes1™	No information
AV300	CaviWipes™ 2.0	No information
AV300	CaviWipes™ HP	No information
AV300	CaviWipes™ Bleach	No information
AV400	CaviCide™	No information
AV400	CaviWipes™	No information
AV400	Envirocide™	No information
AV400	CaviCide1™	Not recommended for use

NEC		
Equipment	Metrex Product	Compatibility Status
NEC Monitors	CaviCide™	No information
NEC Monitors	CaviWipes™	No information
NEC Monitors	Envirocide™	No information
NEC Monitors	CaviCide1™	No information
NEC Monitors	CaviWipes1™	No information
NEC Monitors	CaviWipes™ 2.0	No information
NEC Monitors	CaviWipes™ HP	No information
NEC Monitors	CaviWipes™ Bleach	No information