

# Material Compatibility Status

| ACIST Medical Systems                  |                   |                         |
|----------------------------------------|-------------------|-------------------------|
| Equipment                              | Metrex Product    | Compatibility Status    |
| Angiographic Contrast Injection System | CaviCide™         | Compatible              |
| Angiographic Contrast Injection System | CaviWipes™        | Compatible              |
| Angiographic Contrast Injection System | Envirocide™       | Compatible              |
| Angiographic Contrast Injection System | CaviCide1™        | Not recommended for use |
| Angiographic Contrast Injection System | CaviWipes1™       | Not recommended for use |
| Angiographic Contrast Injection System | CaviWipes™ 2.0    | Compatible              |
| Angiographic Contrast Injection System | CaviWipes™ HP     | Not recommended for use |
| Angiographic Contrast Injection System | CaviWipes™ Bleach | Not recommended for use |

| Physio-Control |                |                      |
|----------------|----------------|----------------------|
| Equipment      | Metrex Product | Compatibility Status |
| LifePak 12     | CaviCide™      | Compatible           |
| LifePak 12     | CaviWipes™     | Compatible           |
| LifePak 12     | Envirocide™    | Compatible           |
| LifePak 12     | CaviCide1™     | Compatible           |
| LifePak 12     | CaviWipes1™    | Compatible           |

| <b>Equipment</b>                          | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------------------------|-----------------------|-----------------------------|
| LifePak 12                                | CaviWipes™ 2.0        | Compatible                  |
| LifePak 12                                | CaviWipes™ HP         | No information              |
| LifePak 12                                | CaviWipes™ Bleach     | No information              |
| LifePak 15                                | CaviCide™             | Compatible                  |
| LifePak 15                                | CaviWipes™            | Compatible                  |
| LifePak 15                                | Envirocide™           | Compatible                  |
| LifePak 15                                | CaviCide1™            | Compatible                  |
| LifePak 15                                | CaviWipes1™           | Compatible                  |
| LifePak 15                                | CaviWipes™ 2.0        | Compatible                  |
| LifePak 15                                | CaviWipes™ HP         | No information              |
| LifePak 15                                | CaviWipes™ Bleach     | No information              |
| LifePak CR Plus Case Only Disposable Pads | CaviCide™             | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes™            | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | Envirocide™           | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviCide1™            | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes1™           | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes™ 2.0        | Compatible                  |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes™ HP         | No information              |
| LifePak CR Plus Case Only Disposable Pads | CaviWipes™ Bleach     | No information              |
| LifePak Express Case only disposable pads | CaviCide™             | Compatible                  |

| <b>Equipment</b>                          | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------------------------|-----------------------|-----------------------------|
| LifePak Express Case only disposable pads | CaviWipes™            | Compatible                  |
| LifePak Express Case only disposable pads | Envirocide™           | Compatible                  |
| LifePak Express Case only disposable pads | CaviCide1™            | Compatible                  |
| LifePak Express Case only disposable pads | CaviWipes1™           | Compatible                  |
| LifePak Express Case only disposable pads | CaviWipes™ 2.0        | Compatible                  |
| LifePak Express Case only disposable pads | CaviWipes™ HP         | No information              |
| LifePak Express Case only disposable pads | CaviWipes™ Bleach     | No information              |