

Material Compatibility Status

Echosens		
Equipment	Metrex Product	Compatibility Status
Fibroscan Mini 430	CaviWipes1™	No information
Fibroscan Mini 430	CaviWipes™ 2.0	No information
Fibroscan Mini 430	CaviWipes™ HP	Not recommended for use
Fibroscan Mini 430	CaviWipes™ Bleach	Not recommended for use
Fibroscan Series	CaviCide™	No information
Fibroscan Series	CaviWipes™	No information
Fibroscan Series	Envirocide™	No information
Fibroscan Series	CaviCide1™	No information
Fibroscan Series	CaviWipes1™	No information
Fibroscan Series	CaviWipes™ 2.0	No information
Fibroscan Series	CaviWipes™ HP	Not recommended for use
Fibroscan Series	CaviWipes™ Bleach	Not recommended for use
New Fibroscan GO	CaviCide™	No information
New Fibroscan GO	CaviWipes™	No information
New Fibroscan GO	Envirocide™	No information
New Fibroscan GO	CaviCide1™	No information
New Fibroscan GO	CaviWipes1™	No information
New Fibroscan GO	CaviWipes™ 2.0	No information
New Fibroscan GO	CaviWipes™ HP	Not recommended for use
New Fibroscan GO	CaviWipes™ Bleach	Not recommended for use
Fibroscan Expert 630	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Fibroscan Expert 630	CaviWipes™	No information
Fibroscan Expert 630	Envirocide™	No information
Fibroscan Expert 630	CaviCide1™	No information
Fibroscan Expert 630	CaviWipes1™	No information
Fibroscan Expert 630	CaviWipes™ 2.0	No information
Fibroscan Expert 630	CaviWipes™ HP	Not recommended for use
Fibroscan Expert 630	CaviWipes™ Bleach	Not recommended for use
Fibroscan compact 530	CaviCide™	No information
Fibroscan compact 530	CaviWipes™	No information
Fibroscan compact 530	Envirocide™	No information
Fibroscan compact 530	CaviCide1™	No information
Fibroscan compact 530	CaviWipes1™	No information
Fibroscan compact 530	CaviWipes™ 2.0	No information
Fibroscan compact 530	CaviWipes™ HP	Not recommended for use
Fibroscan compact 530	CaviWipes™ Bleach	Not recommended for use
Fibroscan Mini+430	CaviCide™	No information
Fibroscan Mini+430	CaviWipes™	No information
Fibroscan Mini+430	Envirocide™	No information
Fibroscan Mini+430	CaviCide1™	No information
Fibroscan Mini+430	CaviWipes1™	No information
Fibroscan Mini+430	CaviWipes™ 2.0	No information
Fibroscan Mini+430	CaviWipes™ HP	Not recommended for use
Fibroscan Mini+430	CaviWipes™ Bleach	Not recommended for use
Fibroscan Mini 430	CaviCide™	No information
Fibroscan Mini 430	CaviWipes™	No information
Fibroscan Mini 430	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Fibrosan Mini 430	CaviCide1™	No information

Fresenius		
Equipment	Metrex Product	Compatibility Status
CRM III Outer Casing	CaviWipes1™	No information
CRM III Outer Casing	CaviWipes™ 2.0	No information
CRM III Outer Casing	CaviWipes™ HP	No information
CRM III Outer Casing	CaviWipes™ Bleach	No information
GranuFlo Concentrate Dissolution Unit II	CaviCide™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™	No information
GranuFlo Concentrate Dissolution Unit II	Envirocide™	No information
GranuFlo Concentrate Dissolution Unit II	CaviCide1™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes1™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ 2.0	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ HP	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ Bleach	Compatible
Liberty Cycler	CaviCide™	Not recommended for use
Liberty Cycler	CaviWipes™	Not recommended for use
Liberty Cycler	Envirocide™	Not recommended for use
Liberty Cycler	CaviCide1™	Not recommended for use
Liberty Cycler	CaviWipes1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Liberty Cyclor	CaviWipes™ 2.0	Not recommended for use
Liberty Cyclor	CaviWipes™ HP	No information
Liberty Cyclor	CaviWipes™ Bleach	Not recommended for use
Newton Cyclor	CaviCide™	Not recommended for use
Newton Cyclor	CaviWipes™	Not recommended for use
Newton Cyclor	Envirocide™	Not recommended for use
Newton Cyclor	CaviCide1™	Not recommended for use
Newton Cyclor	CaviWipes1™	Not recommended for use
Newton Cyclor	CaviWipes™ 2.0	Not recommended for use
Newton Cyclor	CaviWipes™ HP	No information
Newton Cyclor	CaviWipes™ Bleach	Not recommended for use
2008K	CaviCide™	Not recommended for use
2008K	CaviWipes™	Not recommended for use
2008K	Envirocide™	Not recommended for use
2008K	CaviCide1™	Not recommended for use
2008K	CaviWipes1™	Not recommended for use
2008K	CaviWipes™ 2.0	Not recommended for use
2008K	CaviWipes™ HP	No information
2008K	CaviWipes™ Bleach	Not recommended for use
Crit-Line Monitor III Display	CaviCide™	No information
Crit-Line Monitor III Display	CaviWipes™	No information
Crit-Line Monitor III Display	Envirocide™	No information
Crit-Line Monitor III Display	CaviCide1™	No information
Crit-Line Monitor III Display	CaviWipes1™	No information
Crit-Line Monitor III Display	CaviWipes™ 2.0	No information
Crit-Line Monitor III Display	CaviWipes™ HP	No information
Crit-Line Monitor III Display	CaviWipes™ Bleach	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Crit-Line Monitor III Sensor clip	CaviCide™	Compatible
Crit-Line Monitor III Sensor clip	CaviWipes™	Compatible
Crit-Line Monitor III Sensor clip	Envirocide™	Compatible
Crit-Line Monitor III Sensor clip	CaviCide1™	No information
Crit-Line Monitor III Sensor clip	CaviWipes1™	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ 2.0	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ HP	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ Bleach	Compatible
CRM III Outer Casing	CaviCide™	Compatible
CRM III Outer Casing	CaviWipes™	Compatible
CRM III Outer Casing	Envirocide™	Compatible
CRM III Outer Casing	CaviCide1™	No information