

Material Compatibility Status

Hospira		
Equipment	Metrex Product	Compatibility Status
Plum A+3	CaviCide™	No information
Plum A+3	CaviWipes™	No information
Plum A+3	Envirocide™	No information
Plum A+3	CaviCide1™	No information
Plum A+3	CaviWipes1™	No information
Plum A+3	CaviWipes™ 2.0	No information
Plum A+3	CaviWipes™ HP	No information
Plum A+3	CaviWipes™ Bleach	Compatible
Plum XL3	CaviCide™	No information
Plum XL3	CaviWipes™	No information
Plum XL3	Envirocide™	No information
Plum XL3	CaviCide1™	No information
Plum XL3	CaviWipes1™	No information
Plum XL3	CaviWipes™ 2.0	No information
Plum XL3	CaviWipes™ HP	No information
Plum XL3	CaviWipes™ Bleach	Compatible
Plum XL3 Micro/Macro	CaviCide™	No information
Plum XL3 Micro/Macro	CaviWipes™	No information
Plum XL3 Micro/Macro	Envirocide™	No information
Plum XL3 Micro/Macro	CaviCide1™	No information
Plum XL3 Micro/Macro	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Plum XL3 Micro/Macro	CaviWipes™ 2.0	No information
Plum XL3 Micro/Macro	CaviWipes™ HP	No information
Plum XL3 Micro/Macro	CaviWipes™ Bleach	Compatible
Sapphire Pump	CaviCide™	No information
Sapphire Pump	CaviWipes™	No information
Sapphire Pump	Envirocide™	No information
Sapphire Pump	CaviCide1™	No information
Sapphire Pump	CaviWipes1™	No information
Sapphire Pump	CaviWipes™ 2.0	No information
Sapphire Pump	CaviWipes™ HP	No information
Sapphire Pump	CaviWipes™ Bleach	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviCide™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	Envirocide™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes1™	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™ 2.0	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™ HP	No information
Symbiq One Channel Infuser List number16026-04-79 & above, Symbiq Two Channel Infuser List number 16027- 04-79& above	CaviWipes™ Bleach	Compatible
LifeCare PCA List 2709-04	CaviCide™	No information
LifeCare PCA List 2709-04	CaviWipes™	No information
LifeCare PCA List 2709-04	Envirocide™	No information
LifeCare PCA List 2709-04	CaviCide1™	No information
LifeCare PCA List 2709-04	CaviWipes1™	No information
LifeCare PCA List 2709-04	CaviWipes™ 2.0	No information
LifeCare PCA List 2709-04	CaviWipes™ HP	No information
LifeCare PCA List 2709-04	CaviWipes™ Bleach	Compatible
Omni Flow 4000Plus	CaviCide™	No information
Omni Flow 4000Plus	CaviWipes™	No information
Omni Flow 4000Plus	Envirocide™	No information
Omni Flow 4000Plus	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Omni Flow 4000Plus	CaviWipes1™	No information
Omni Flow 4000Plus	CaviWipes™ 2.0	No information
Omni Flow 4000Plus	CaviWipes™ HP	No information
Omni Flow 4000Plus	CaviWipes™ Bleach	Compatible
Plum 360	CaviCide™	No information
Plum 360	CaviWipes™	No information
Plum 360	Envirocide™	No information
Plum 360	CaviCide1™	No information
Plum 360	CaviWipes1™	No information
Plum 360	CaviWipes™ 2.0	No information
Plum 360	CaviWipes™ HP	No information
Plum 360	CaviWipes™ Bleach	Compatible
Plum A+	CaviCide™	No information
Plum A+	CaviWipes™	No information
Plum A+	Envirocide™	No information
Plum A+	CaviCide1™	No information
Plum A+	CaviWipes1™	No information
Plum A+	CaviWipes™ 2.0	No information
Plum A+	CaviWipes™ HP	No information
Plum A+	CaviWipes™ Bleach	Compatible
Plum A+ W/Hospira MedNet	CaviCide™	No information
Plum A+ W/Hospira MedNet	CaviWipes™	No information
Plum A+ W/Hospira MedNet	Envirocide™	No information
Plum A+ W/Hospira MedNet	CaviCide1™	No information
Plum A+ W/Hospira MedNet	CaviWipes1™	No information
Plum A+ W/Hospira MedNet	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Plum A+ W/Hospira MedNet	CaviWipes™ HP	No information
Plum A+ W/Hospira MedNet	CaviWipes™ Bleach	Compatible
Acclaim Encore	CaviCide™	No information
Acclaim Encore	CaviWipes™	No information
Acclaim Encore	Envirocide™	No information
Acclaim Encore	CaviCide1™	No information
Acclaim Encore	CaviWipes1™	No information
Acclaim Encore	CaviWipes™ 2.0	No information
Acclaim Encore	CaviWipes™ HP	No information
Acclaim Encore	CaviWipes™ Bleach	Compatible
GemStar	CaviCide™	No information
GemStar	CaviWipes™	No information
GemStar	Envirocide™	No information
GemStar	CaviCide1™	No information
GemStar	CaviWipes1™	No information
GemStar	CaviWipes™ 2.0	No information
GemStar	CaviWipes™ HP	No information
GemStar	CaviWipes™ Bleach	Compatible

RMS Medical Products		
Equipment	Metrex Product	Compatibility Status
Freedom 60 Syringe Infusion system	CaviCide™	Compatible
Freedom 60 Syringe Infusion system	CaviWipes™	Compatible
Freedom 60 Syringe Infusion system	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Freedom 60 Syringe Infusion system	CaviCide1™	No information
Freedom 60 Syringe Infusion system	CaviWipes1™	No information
Freedom 60 Syringe Infusion system	CaviWipes™ 2.0	No information
Freedom 60 Syringe Infusion system	CaviWipes™ HP	No information
Freedom 60 Syringe Infusion system	CaviWipes™ Bleach	Compatible