

# Material Compatibility Status

| LifeScan               |                   |                         |
|------------------------|-------------------|-------------------------|
| Equipment              | Metrex Product    | Compatibility Status    |
| One Touch Delica       | CaviCide™         | No information          |
| One Touch Delica       | CaviWipes™        | No information          |
| One Touch Delica       | Envirocide™       | No information          |
| One Touch Delica       | CaviCide1™        | No information          |
| One Touch Delica       | CaviWipes1™       | No information          |
| One Touch Delica       | CaviWipes™ 2.0    | No information          |
| One Touch Delica       | CaviWipes™ HP     | No information          |
| One Touch Delica       | CaviWipes™ Bleach | Compatible              |
| Onetouch Ping          | CaviCide™         | No information          |
| Onetouch Ping          | CaviWipes™        | No information          |
| Onetouch Ping          | Envirocide™       | No information          |
| Onetouch Ping          | CaviCide1™        | No information          |
| Onetouch Ping          | CaviWipes1™       | No information          |
| Onetouch Ping          | CaviWipes™ 2.0    | No information          |
| Onetouch Ping          | CaviWipes™ HP     | No information          |
| Onetouch Ping          | CaviWipes™ Bleach | Not recommended for use |
| Onetouch Verio systems | CaviWipes1™       | No information          |
| Onetouch Verio systems | CaviWipes™ 2.0    | No information          |
| Onetouch Verio systems | CaviWipes™ HP     | No information          |
| Onetouch Verio systems | CaviWipes™ Bleach | No information          |
| Surestep Flex          | CaviCide™         | No information          |

| <b>Equipment</b>   | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|--------------------|-----------------------|-----------------------------|
| Surestep Flex      | CaviWipes™            | No information              |
| Surestep Flex      | Envirocide™           | No information              |
| Surestep Flex      | CaviCide1™            | No information              |
| Surestep Flex      | CaviWipes1™           | No information              |
| Surestep Flex      | CaviWipes™ 2.0        | No information              |
| Surestep Flex      | CaviWipes™ HP         | No information              |
| Surestep Flex      | CaviWipes™ Bleach     | Compatible                  |
| Surestep Meters    | CaviCide™             | No information              |
| Surestep Meters    | CaviWipes™            | No information              |
| Surestep Meters    | Envirocide™           | No information              |
| Surestep Meters    | CaviCide1™            | No information              |
| Surestep Meters    | CaviWipes1™           | No information              |
| Surestep Meters    | CaviWipes™ 2.0        | No information              |
| Surestep Meters    | CaviWipes™ HP         | No information              |
| Surestep Meters    | CaviWipes™ Bleach     | Compatible                  |
| Onetouch Ultra 2   | CaviCide™             | No information              |
| Onetouch Ultra 2   | CaviWipes™            | No information              |
| Onetouch Ultra 2   | Envirocide™           | No information              |
| Onetouch Ultra 2   | CaviCide1™            | No information              |
| Onetouch Ultra 2   | CaviWipes1™           | No information              |
| Onetouch Ultra 2   | CaviWipes™ 2.0        | No information              |
| Onetouch Ultra 2   | CaviWipes™ HP         | No information              |
| Onetouch Ultra 2   | CaviWipes™ Bleach     | Compatible                  |
| Onetouch Ultramini | CaviCide™             | No information              |
| Onetouch Ultramini | CaviWipes™            | No information              |
| Onetouch Ultramini | Envirocide™           | No information              |

| <b>Equipment</b>       | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------------|-----------------------|-----------------------------|
| Onetouch Ultramini     | CaviCide1™            | No information              |
| Onetouch Ultramini     | CaviWipes1™           | No information              |
| Onetouch Ultramini     | CaviWipes™ 2.0        | No information              |
| Onetouch Ultramini     | CaviWipes™ HP         | No information              |
| Onetouch Ultramini     | CaviWipes™ Bleach     | Compatible                  |
| Onetouch Verio systems | CaviCide™             | No information              |
| Onetouch Verio systems | CaviWipes™            | No information              |
| Onetouch Verio systems | Envirocide™           | No information              |
| Onetouch Verio systems | CaviCide1™            | No information              |

| <b>Panasonic</b> |                       |                             |
|------------------|-----------------------|-----------------------------|
| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |
| ToughBook        | CaviCide™             | Compatible                  |
| ToughBook        | CaviWipes™            | Compatible                  |
| ToughBook        | Envirocide™           | Compatible                  |
| ToughBook        | CaviCide1™            | No information              |
| ToughBook        | CaviWipes1™           | No information              |
| ToughBook        | CaviWipes™ 2.0        | No information              |
| ToughBook        | CaviWipes™ HP         | No information              |
| ToughBook        | CaviWipes™ Bleach     | No information              |
| ToughBook        | CaviCide™             | Compatible                  |
| ToughBook        | CaviWipes™            | Compatible                  |
| ToughBook        | Envirocide™           | Compatible                  |
| ToughBook        | CaviCide1™            | No information              |
| ToughBook        | CaviWipes1™           | No information              |
| ToughBook        | CaviWipes™ 2.0        | No information              |
| ToughBook        | CaviWipes™ HP         | No information              |

| <b>Equipment</b>     | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|----------------------|-----------------------|-----------------------------|
| ToughBook            | CaviWipes™ Bleach     | No information              |
| Upper arm BP Monitor | CaviCide™             | No information              |
| Upper arm BP Monitor | CaviWipes™            | No information              |
| Upper arm BP Monitor | Envirocide™           | No information              |
| Upper arm BP Monitor | CaviCide1™            | No information              |
| Upper arm BP Monitor | CaviWipes1™           | No information              |
| Upper arm BP Monitor | CaviWipes™ 2.0        | No information              |
| Upper arm BP Monitor | CaviWipes™ HP         | No information              |
| Upper arm BP Monitor | CaviWipes™ Bleach     | No information              |
| Wrist BP Monitor     | CaviCide™             | No information              |
| Wrist BP Monitor     | CaviWipes™            | No information              |
| Wrist BP Monitor     | Envirocide™           | No information              |
| Wrist BP Monitor     | CaviCide1™            | No information              |
| Wrist BP Monitor     | CaviWipes1™           | No information              |
| Wrist BP Monitor     | CaviWipes™ 2.0        | No information              |
| Wrist BP Monitor     | CaviWipes™ HP         | No information              |
| Wrist BP Monitor     | CaviWipes™ Bleach     | No information              |