

Material Compatibility Status

LifeScan		
Equipment	Metrex Product	Compatibility Status
One Touch Delica	CaviCide™	No information
One Touch Delica	CaviWipes™	No information
One Touch Delica	Envirocide™	No information
One Touch Delica	CaviCide1™	No information
One Touch Delica	CaviWipes1™	No information
One Touch Delica	CaviWipes™ 2.0	No information
One Touch Delica	CaviWipes™ HP	No information
One Touch Delica	CaviWipes™ Bleach	Compatible
Onetouch Ping	CaviCide™	No information
Onetouch Ping	CaviWipes™	No information
Onetouch Ping	Envirocide™	No information
Onetouch Ping	CaviCide1™	No information
Onetouch Ping	CaviWipes1™	No information
Onetouch Ping	CaviWipes™ 2.0	No information
Onetouch Ping	CaviWipes™ HP	No information
Onetouch Ping	CaviWipes™ Bleach	Not recommended for use
Onetouch Verio systems	CaviWipes1™	No information
Onetouch Verio systems	CaviWipes™ 2.0	No information
Onetouch Verio systems	CaviWipes™ HP	No information
Onetouch Verio systems	CaviWipes™ Bleach	No information
Surestep Flex	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Surestep Flex	CaviWipes™	No information
Surestep Flex	Envirocide™	No information
Surestep Flex	CaviCide1™	No information
Surestep Flex	CaviWipes1™	No information
Surestep Flex	CaviWipes™ 2.0	No information
Surestep Flex	CaviWipes™ HP	No information
Surestep Flex	CaviWipes™ Bleach	Compatible
Surestep Meters	CaviCide™	No information
Surestep Meters	CaviWipes™	No information
Surestep Meters	Envirocide™	No information
Surestep Meters	CaviCide1™	No information
Surestep Meters	CaviWipes1™	No information
Surestep Meters	CaviWipes™ 2.0	No information
Surestep Meters	CaviWipes™ HP	No information
Surestep Meters	CaviWipes™ Bleach	Compatible
Onetouch Ultra 2	CaviCide™	No information
Onetouch Ultra 2	CaviWipes™	No information
Onetouch Ultra 2	Envirocide™	No information
Onetouch Ultra 2	CaviCide1™	No information
Onetouch Ultra 2	CaviWipes1™	No information
Onetouch Ultra 2	CaviWipes™ 2.0	No information
Onetouch Ultra 2	CaviWipes™ HP	No information
Onetouch Ultra 2	CaviWipes™ Bleach	Compatible
Onetouch Ultramini	CaviCide™	No information
Onetouch Ultramini	CaviWipes™	No information
Onetouch Ultramini	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Onetouch Ultramini	CaviCide1™	No information
Onetouch Ultramini	CaviWipes1™	No information
Onetouch Ultramini	CaviWipes™ 2.0	No information
Onetouch Ultramini	CaviWipes™ HP	No information
Onetouch Ultramini	CaviWipes™ Bleach	Compatible
Onetouch Verio systems	CaviCide™	No information
Onetouch Verio systems	CaviWipes™	No information
Onetouch Verio systems	Envirocide™	No information
Onetouch Verio systems	CaviCide1™	No information

RMS Medical Products		
Equipment	Metrex Product	Compatibility Status
Freedom 60 Syringe Infusion system	CaviCide™	Compatible
Freedom 60 Syringe Infusion system	CaviWipes™	Compatible
Freedom 60 Syringe Infusion system	Envirocide™	Compatible
Freedom 60 Syringe Infusion system	CaviCide1™	No information
Freedom 60 Syringe Infusion system	CaviWipes1™	No information
Freedom 60 Syringe Infusion system	CaviWipes™ 2.0	No information
Freedom 60 Syringe Infusion system	CaviWipes™ HP	No information
Freedom 60 Syringe Infusion system	CaviWipes™ Bleach	Compatible