

Material Compatibility Status

MediStim		
Equipment	Metrex Product	Compatibility Status
PV	CaviCide™	Compatible
PV	CaviWipes™	Compatible
PV	Envirocide™	Compatible
PV	CaviCide1™	No information
PV	CaviWipes1™	No information
PV	CaviWipes™ 2.0	No information
PV	CaviWipes™ HP	No information
PV	CaviWipes™ Bleach	No information
PD	CaviCide™	Compatible
PD	CaviWipes™	Compatible
PD	Envirocide™	Compatible
PD	CaviCide1™	No information
PD	CaviWipes1™	No information
PD	CaviWipes™ 2.0	No information
PD	CaviWipes™ HP	No information
PD	CaviWipes™ Bleach	No information