

Material Compatibility Status

Moldex Hearing & Respiratory		
Equipment	Metrex Product	Compatibility Status
7000 Series	CaviCide™	No information
7000 Series	CaviWipes™	No information
7000 Series	Envirocide™	No information
7000 Series	CaviCide1™	No information
7000 Series	CaviWipes1™	No information
7000 Series	CaviWipes™ 2.0	No information
7000 Series	CaviWipes™ HP	No information
7000 Series	CaviWipes™ Bleach	No information

RMS Medical Products		
Equipment	Metrex Product	Compatibility Status
Freedom 60 Syringe Infusion system	CaviCide™	Compatible
Freedom 60 Syringe Infusion system	CaviWipes™	Compatible
Freedom 60 Syringe Infusion system	Envirocide™	Compatible
Freedom 60 Syringe Infusion system	CaviCide1™	No information
Freedom 60 Syringe Infusion system	CaviWipes1™	No information
Freedom 60 Syringe Infusion system	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Freedom 60 Syringe Infusion system	CaviWipes™ HP	No information
Freedom 60 Syringe Infusion system	CaviWipes™ Bleach	Compatible