

Material Compatibility Status

| Spacelabs | | |
|------------------------------------|-------------------|----------------------|
| Equipment | Metrex Product | Compatibility Status |
| Tri-Cuff | CaviWipes™ Bleach | No information |
| TRU CUFF | CaviWipes™ Bleach | Compatible |
| TRU LINK | CaviWipes™ Bleach | Compatible |
| Ultaview SL | CaviWipes™ Bleach | Compatible |
| Ultraview Care Network | CaviWipes™ Bleach | No information |
| Ultraview Care Network/Monitor | CaviWipes™ Bleach | No information |
| Ultraview Care Network/Transducers | CaviWipes™ Bleach | No information |
| Ultraview Care Network/Sensors | CaviWipes™ Bleach | Compatible |