

Material Compatibility Status

Resmed		
Equipment	Metrex Product	Compatibility Status
Astral 100	CaviWipes1™	No information
Astral 100	CaviWipes™ 2.0	No information
Astral 100	CaviWipes™ HP	No information
Astral 100	CaviWipes™ Bleach	Compatible
Astral 150	CaviCide™	Compatible
Astral 150	CaviWipes™	Compatible
Astral 150	Envirocide™	Compatible
Astral 150	CaviCide1™	No information
Astral 150	CaviWipes1™	No information
Astral 150	CaviWipes™ 2.0	No information
Astral 150	CaviWipes™ HP	No information
Astral 150	CaviWipes™ Bleach	Compatible
Astral 100	CaviCide™	Compatible
Astral 100	CaviWipes™	Compatible
Astral 100	Envirocide™	Compatible
Astral 100	CaviCide1™	No information