

Material Compatibility Status

RMS Medical Products		
Equipment	Metrex Product	Compatibility Status
Freedom 60 Syringe Infusion system	CaviCide™	Compatible
Freedom 60 Syringe Infusion system	CaviWipes™	Compatible
Freedom 60 Syringe Infusion system	Envirocide™	Compatible
Freedom 60 Syringe Infusion system	CaviCide1™	No information
Freedom 60 Syringe Infusion system	CaviWipes1™	No information
Freedom 60 Syringe Infusion system	CaviWipes™ 2.0	No information
Freedom 60 Syringe Infusion system	CaviWipes™ HP	No information
Freedom 60 Syringe Infusion system	CaviWipes™ Bleach	Compatible