

Material Compatibility Status

Dentlight		
Equipment	Metrex Product	Compatibility Status
Fusion 5 Handpiece & Charger Base	CaviCide™	Compatible
Fusion 5 Handpiece & Charger Base	CaviWipes™	Compatible
Fusion 5 Handpiece & Charger Base	Envirocide™	Compatible
Fusion 5 Handpiece & Charger Base	CaviCide1™	No information
Fusion 5 Handpiece & Charger Base	CaviWipes1™	No information
Fusion 5 Handpiece & Charger Base	CaviWipes™ 2.0	No information
Fusion 5 Handpiece & Charger Base	CaviWipes™ HP	No information
Fusion 5 Handpiece & Charger Base	CaviWipes™ Bleach	No information