

# Material Compatibility Status

3M		
Equipment	Metrex Product	Compatibility Status
M-Series Headgear	CaviCide™	Not recommended for use
M-Series Headgear	CaviWipes™	Not recommended for use
M-Series Headgear	Envirocide™	Not recommended for use
M-Series Headgear	CaviCide1™	Not recommended for use
M-Series Headgear	CaviWipes1™	Not recommended for use
M-Series Headgear	CaviWipes™ 2.0	Not recommended for use

A-dec		
Equipment	Metrex Product	Compatibility Status
Preference Collection Model 5562 Laminate & quartz	CaviWipes™ Bleach	Not recommended for use
Preference Collection Model 5580 Laminate & quartz	CaviWipes™ Bleach	Not recommended for use
Inspire Collection: Infills	CaviWipes™ Bleach	Not recommended for use
Inspire Collection: Laminate surfaces	CaviWipes™ Bleach	Not recommended for use
Inspire Collection: Quartz surfaces	CaviWipes™ Bleach	Not recommended for use
Inspire Collection: Solid surfaces	CaviWipes™ Bleach	Not recommended for use
A-dec 400	CaviWipes™ Bleach	Not recommended for use
A-dec 500	CaviWipes™ Bleach	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
A-dec 511b	CaviWipes™ HP	Not recommended for use
A-dec 511b	CaviWipes™ Bleach	Not recommended for use
A-dec 300	CaviWipes™ Bleach	Not recommended for use
Preference Collection Model 5562 Treatment console insert / solid surface	CaviWipes™ Bleach	Not recommended for use
Preference Collection Model 5580 Treatment console / solid surface	CaviWipes™ Bleach	Not recommended for use
Inspire 391 Treatment Console	CaviWipes™ Bleach	Not recommended for use
Inspire 591 A Treatment Console	CaviWipes™ Bleach	Not recommended for use
Inspire 592 A Central Console	CaviWipes™ Bleach	Not recommended for use
A-dec 311	CaviWipes™ HP	Not recommended for use
A-dec 311	CaviWipes™ Bleach	Not recommended for use
A-dec 411	CaviWipes™ HP	Not recommended for use
A-dec 411	CaviWipes™ Bleach	Not recommended for use
Inspire 593 a Side Console	CaviWipes™ Bleach	Not recommended for use
Inspire 594 A Sterilization Center	CaviWipes™ Bleach	Not recommended for use
Inspire 595 A Wall Mounted Cabinet/Xray Storage Insert	CaviWipes™ Bleach	Not recommended for use
A-dec 511	CaviWipes™ HP	Not recommended for use
A-dec 511	CaviWipes™ Bleach	Not recommended for use
A-dec Upholstery	CaviWipes™ HP	Not recommended for use
A-dec Upholstery	CaviWipes™ Bleach	Not recommended for use
Model 500 LED	CaviWipes™ Bleach	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Model 300 LED	CaviWipes™ Bleach	Not recommended for use
Model 421	CaviWipes™ HP	Not recommended for use
Model 421	CaviWipes™ Bleach	Not recommended for use
Model 422	CaviWipes™ HP	Not recommended for use
Model 422	CaviWipes™ Bleach	Not recommended for use
Model 521	CaviWipes™ HP	Not recommended for use
Model 521	CaviWipes™ Bleach	Not recommended for use
Model 522	CaviWipes™ HP	Not recommended for use
Model 522	CaviWipes™ Bleach	Not recommended for use

<b>Abbott Healthcare</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Precision Xceed Pro	CaviCide1™	Not recommended for use
Precision Xceed Pro	CaviWipes1™	Not recommended for use
Precision Xceed Pro	CaviWipes™ HP	Not recommended for use
Precision Xceed Pro	CaviWipes™ Bleach	Not recommended for use
Precision Xtra	CaviCide™	Not recommended for use
Precision Xtra	CaviWipes™	Not recommended for use
Precision Xtra	Envirocide™	Not recommended for use
Precision Xtra	CaviCide1™	Not recommended for use
Precision Xtra	CaviWipes1™	Not recommended for use
Precision Xtra	CaviWipes™ 2.0	Not recommended for use
Precision Xtra	CaviWipes™ HP	Not recommended for use
Flexible Quantum	CaviCide™	Not recommended for use
Flexible Quantum	CaviWipes™	Not recommended for use
Flexible Quantum	Envirocide™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Flexible Quantum	CaviCide1™	Not recommended for use
Flexible Quantum	CaviWipes1™	Not recommended for use
Flexible Quantum	CaviWipes™ 2.0	Not recommended for use
Flexible Quantum	CaviWipes™ HP	Not recommended for use
Freestyle Flash	CaviCide™	Not recommended for use
Freestyle Flash	CaviWipes™	Not recommended for use
Freestyle Flash	Envirocide™	Not recommended for use
Freestyle Flash	CaviCide1™	Not recommended for use
Freestyle Flash	CaviWipes1™	Not recommended for use
Freestyle Flash	CaviWipes™ 2.0	Not recommended for use
Freestyle Flash	CaviWipes™ HP	Not recommended for use
Freestyle Freedom	CaviCide™	Not recommended for use
Freestyle Freedom	CaviWipes™	Not recommended for use
Freestyle Freedom	Envirocide™	Not recommended for use
Freestyle Freedom	CaviCide1™	Not recommended for use
Freestyle Freedom	CaviWipes1™	Not recommended for use
Freestyle Freedom	CaviWipes™ 2.0	Not recommended for use
Freestyle Freedom	CaviWipes™ HP	Not recommended for use
Freestyle Lite	CaviWipes1™	Not recommended for use
Freestyle Lite	CaviWipes™ 2.0	Not recommended for use
Freestyle Lite	CaviWipes™ HP	Not recommended for use
Optimum Xceed	CaviCide™	Not recommended for use
Optimum Xceed	CaviWipes™	Not recommended for use
Optimum Xceed	Envirocide™	Not recommended for use
Optimum Xceed	CaviCide1™	Not recommended for use
Optimum Xceed	CaviWipes1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Optimum Xceed	CaviWipes™ 2.0	Not recommended for use
Optimum Xceed	CaviWipes™ HP	Not recommended for use
Patrol ePump	CaviCide™	Not recommended for use
Patrol ePump	CaviWipes™	Not recommended for use
Patrol ePump	Envirocide™	Not recommended for use
Patrol ePump	CaviCide1™	Not recommended for use
Patrol ePump	CaviWipes1™	Not recommended for use
Patrol ePump	CaviWipes™ 2.0	Not recommended for use
Patrol ePump	CaviWipes™ HP	Not recommended for use
FreeStyle Freedom Lite	CaviCide™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™	Not recommended for use
FreeStyle Freedom Lite	Envirocide™	Not recommended for use
FreeStyle Freedom Lite	CaviCide1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes1™	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ 2.0	Not recommended for use
FreeStyle Freedom Lite	CaviWipes™ HP	Not recommended for use
Freestyle Insulinx	CaviCide™	Not recommended for use
Freestyle Insulinx	CaviWipes™	Not recommended for use
Freestyle Insulinx	Envirocide™	Not recommended for use
Freestyle Insulinx	CaviCide1™	Not recommended for use
Freestyle Insulinx	CaviWipes1™	Not recommended for use
Freestyle Insulinx	CaviWipes™ 2.0	Not recommended for use
Freestyle Insulinx	CaviWipes™ HP	Not recommended for use
Freestyle Insulinx	CaviWipes™ Bleach	Not recommended for use
Freestyle Lite	CaviCide™	Not recommended for use
Freestyle Lite	CaviWipes™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Freestyle Lite	Envirocide™	Not recommended for use
Freestyle Lite	CaviCide1™	Not recommended for use

<b>Accuvein</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
AV400	CaviWipes1™	Not recommended for use
AV 300	CaviCide1™	Not recommended for use
AV 300	CaviWipes1™	Not recommended for use
AV400	CaviCide1™	Not recommended for use

<b>ACIST Medical Systems</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Angiographic Contrast Injection System	CaviCide1™	Not recommended for use
Angiographic Contrast Injection System	CaviWipes1™	Not recommended for use
Angiographic Contrast Injection System	CaviWipes™ HP	Not recommended for use
Angiographic Contrast Injection System	CaviWipes™ Bleach	Not recommended for use

<b>AgaMatrix</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Jazz	CaviCide™	Not recommended for use
Jazz	CaviWipes™	Not recommended for use
Jazz	Envirocide™	Not recommended for use
Jazz	CaviCide1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Jazz	CaviWipes1™	Not recommended for use
Jazz	CaviWipes™ 2.0	Not recommended for use
Jazz	CaviWipes™ HP	Not recommended for use
Jazz	CaviWipes™ Bleach	Not recommended for use

<b>AGFA</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
DR 800	CaviWipes™ HP	Not recommended for use
DR 800	CaviWipes™ Bleach	Not recommended for use

<b>Air Techniques</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Provecta S-Pan Ceph imaging	CaviCide™	Not recommended for use
Provecta S-Pan Ceph imaging	CaviWipes™	Not recommended for use
Provecta S-Pan Ceph imaging	Envirocide™	Not recommended for use
Provecta S-Pan Ceph imaging	CaviCide1™	Not recommended for use
Provecta S-Pan Ceph imaging	CaviWipes1™	Not recommended for use

<b>Arizant Healthcare</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Bair Hugger 700 Warming Units	CaviCide1™	Not recommended for use
Bair Hugger 700 Warming Units	CaviWipes1™	Not recommended for use
Bair Hugger 800 Warming Units	CaviCide1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Bair Hugger 800 Warming Units	CaviWipes1™	Not recommended for use
Bair Paws Model 850	CaviWipes™ Bleach	Not recommended for use
Ranger Model 245	CaviCide1™	Not recommended for use
Ranger Model 245	CaviWipes1™	Not recommended for use
Ranger Model 247	CaviCide1™	Not recommended for use
Ranger Model 247	CaviWipes1™	Not recommended for use

<b>Belmont</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
D Series/Model: D-4R	CaviWipes™ HP	Not recommended for use
D Series/Model: D-4R	CaviWipes™ Bleach	Not recommended for use
D Series/Model: D-4RH	CaviWipes™ HP	Not recommended for use
D Series/Model: D-4RH	CaviWipes™ Bleach	Not recommended for use
Bel-20+	CaviWipes™ Bleach	Not recommended for use
Bel-50N	CaviWipes™ Bleach	Not recommended for use
X-Calibur V B50N	CaviWipes™ Bleach	Not recommended for use
Synthetic leather upholstery	CaviWipes™ Bleach	Not recommended for use
Quolis 5500	CaviWipes™ Bleach	Not recommended for use
Bel-Ray II AC	CaviWipes™ Bleach	Not recommended for use
X-Calibur Bel-2N	CaviWipes™ Bleach	Not recommended for use
PROIII 047	CaviWipes™ Bleach	Not recommended for use
X-Caliber BDS	CaviWipes™ Bleach	Not recommended for use
Quolis 5000 Series	CaviWipes™ Bleach	Not recommended for use

<b>Brewer</b>
---------------



<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Halogen 35 Series Model 18100	CaviWipes™ HP	Not recommended for use
Incandescent Series Model 17100	CaviWipes™ HP	Not recommended for use
Model 11220	CaviWipes™ HP	Not recommended for use
Model 11200-1	CaviWipes™ HP	Not recommended for use
Model 11200-2	CaviWipes™ HP	Not recommended for use
All Purpose Stainless Steel	CaviWipes™ Bleach	Not recommended for use
Model 31200	CaviWipes™ HP	Not recommended for use
Model 31220	CaviWipes™ HP	Not recommended for use
Model 21200	CaviWipes™ HP	Not recommended for use

<b>DCI Edge</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Orthodontic chair	CaviWipes™ HP	Not recommended for use
Orthodontic chair	CaviWipes™ Bleach	Not recommended for use

<b>Dentsply</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
SmartLite Focus	CaviWipes1™	Not recommended for use
SmartLite Focus	CaviCide™	Not recommended for use
SmartLite Focus	CaviWipes™	Not recommended for use
SmartLite Focus	Envirocide™	Not recommended for use
SmartLite Focus	CaviCide1™	Not recommended for use

<b>Draeger</b>		
----------------	--	--

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Fabius MRI all respiratory components	CaviWipes™ HP	Not recommended for use

<b>Echosens</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Fibroscan Mini 430	CaviWipes™ HP	Not recommended for use
Fibroscan Mini 430	CaviWipes™ Bleach	Not recommended for use
Fibroscan Series	CaviWipes™ HP	Not recommended for use
Fibroscan Series	CaviWipes™ Bleach	Not recommended for use
New Fibroscan GO	CaviWipes™ HP	Not recommended for use
New Fibroscan GO	CaviWipes™ Bleach	Not recommended for use
Fibroscan Expert 630	CaviWipes™ HP	Not recommended for use
Fibroscan Expert 630	CaviWipes™ Bleach	Not recommended for use
Fibroscan compact 530	CaviWipes™ HP	Not recommended for use
Fibroscan compact 530	CaviWipes™ Bleach	Not recommended for use
Fibroscan Mini+430	CaviWipes™ HP	Not recommended for use
Fibroscan Mini+430	CaviWipes™ Bleach	Not recommended for use

<b>Engle Dental</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Engle 300	CaviWipes™ Bleach	Not recommended for use

<b>Fresenius</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Liberty Cyclor	CaviCide™	Not recommended for use
Liberty Cyclor	CaviWipes™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Liberty Cyclor	Envirocide™	Not recommended for use
Liberty Cyclor	CaviCide1™	Not recommended for use
Liberty Cyclor	CaviWipes1™	Not recommended for use
Liberty Cyclor	CaviWipes™ 2.0	Not recommended for use
Liberty Cyclor	CaviWipes™ Bleach	Not recommended for use
Newton Cyclor	CaviCide™	Not recommended for use
Newton Cyclor	CaviWipes™	Not recommended for use
Newton Cyclor	Envirocide™	Not recommended for use
Newton Cyclor	CaviCide1™	Not recommended for use
Newton Cyclor	CaviWipes1™	Not recommended for use
Newton Cyclor	CaviWipes™ 2.0	Not recommended for use
Newton Cyclor	CaviWipes™ Bleach	Not recommended for use
2008K	CaviCide™	Not recommended for use
2008K	CaviWipes™	Not recommended for use
2008K	Envirocide™	Not recommended for use
2008K	CaviCide1™	Not recommended for use
2008K	CaviWipes1™	Not recommended for use
2008K	CaviWipes™ 2.0	Not recommended for use
2008K	CaviWipes™ Bleach	Not recommended for use
Crit-Line Monitor III Display	CaviWipes™ Bleach	Not recommended for use

<b>FujiFilm</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
C11x	CaviCide1™	Not recommended for use
C11x	CaviWipes1™	Not recommended for use
C35X	CaviCide1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
L38x	CaviWipes1™	Not recommended for use
L38xi	CaviCide1™	Not recommended for use
L38xi	CaviWipes1™	Not recommended for use
Power supply	CaviCide1™	Not recommended for use
Power supply	CaviWipes1™	Not recommended for use
SLAx	CaviCide1™	Not recommended for use
SLAx	CaviWipes1™	Not recommended for use
NanoMaxx Field Case	CaviCide1™	Not recommended for use
NanoMaxx Field Case	CaviWipes1™	Not recommended for use
H-Universal Stand	CaviCide1™	Not recommended for use
C35X	CaviWipes1™	Not recommended for use
rC60xi	CaviCide1™	Not recommended for use
rC60xi	CaviWipes1™	Not recommended for use
ICTx Revision 80 and below	CaviCide1™	Not recommended for use
ICTx Revision 80 and below	CaviWipes1™	Not recommended for use
ICTx Revision 81	CaviCide1™	Not recommended for use
ICTx Revision 81	CaviWipes1™	Not recommended for use
C10-3	CaviCide1™	Not recommended for use
C10-3	CaviWipes1™	Not recommended for use
C11n	CaviCide1™	Not recommended for use
Edge II Stand	CaviWipes1™	Not recommended for use
EDGE II	CaviCide1™	Not recommended for use
EDGE II	CaviWipes1™	Not recommended for use
TTC - Edge	CaviCide1™	Not recommended for use
TTC - Edge	CaviWipes1™	Not recommended for use
HFL50X	CaviCide1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
HFL50X	CaviWipes1™	Not recommended for use
L25x	CaviCide1™	Not recommended for use
L25x	CaviWipes1™	Not recommended for use
L38x	CaviCide1™	Not recommended for use
C11n	CaviWipes1™	Not recommended for use
C11xp	CaviCide1™	Not recommended for use
C11xp	CaviWipes1™	Not recommended for use
C35xp	CaviCide1™	Not recommended for use
C35xp	CaviWipes1™	Not recommended for use
C38xi	CaviCide1™	Not recommended for use
C38xi	CaviWipes1™	Not recommended for use
C5-1	CaviCide1™	Not recommended for use
C5-1	CaviWipes1™	Not recommended for use
C60n	CaviCide1™	Not recommended for use
C60xf	CaviCide1™	Not recommended for use
C60xf	CaviWipes1™	Not recommended for use
C60xi	CaviCide1™	Not recommended for use
C60xi	CaviWipes1™	Not recommended for use
C60xp	CaviCide1™	Not recommended for use
C60xp	CaviWipes1™	Not recommended for use
C8x	CaviCide1™	Not recommended for use
C8x	CaviWipes1™	Not recommended for use
D2xp	CaviCide™	Not recommended for use
D2xp	CaviWipes™	Not recommended for use
D2xp	Envirocide™	Not recommended for use
D2xp	CaviCide1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
D2xp	CaviWipes1™	Not recommended for use
H-Universal Stand	CaviWipes1™	Not recommended for use
HFL38xp	CaviCide1™	Not recommended for use
HFL38xp	CaviWipes1™	Not recommended for use
HFL38x	CaviCide1™	Not recommended for use
HFL38x	CaviWipes1™	Not recommended for use
Dock	CaviCide1™	Not recommended for use
Dock	CaviWipes1™	Not recommended for use
ECG Cable	CaviCide1™	Not recommended for use
ECG Cable	CaviWipes1™	Not recommended for use
Edge II Stand	CaviCide1™	Not recommended for use
HFL38xi	CaviCide1™	Not recommended for use
HFL38xi	CaviWipes1™	Not recommended for use
HFL50xp	CaviCide1™	Not recommended for use
HFL50xp	CaviWipes1™	Not recommended for use
HSL25x	CaviCide1™	Not recommended for use
C60n	CaviWipes1™	Not recommended for use
C60v	CaviCide1™	Not recommended for use
C60v	CaviWipes1™	Not recommended for use
C60x	CaviCide1™	Not recommended for use
C60x	CaviWipes1™	Not recommended for use
IC10-3	CaviCide1™	Not recommended for use
IC10-3	CaviWipes1™	Not recommended for use
ICTn	CaviCide1™	Not recommended for use
ICTn	CaviWipes1™	Not recommended for use
ICTxp Revision 21 and below	CaviCide1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
ICTxp Revision 21 and below	CaviWipes1™	Not recommended for use
ICTxp Revision 22	CaviCide1™	Not recommended for use
ICTxp Revision 22	CaviWipes1™	Not recommended for use
L38V	CaviCide1™	Not recommended for use
L38V	CaviWipes1™	Not recommended for use
L25n	CaviWipes1™	Not recommended for use
L25v	CaviCide1™	Not recommended for use
L25v	CaviWipes1™	Not recommended for use
L25x	CaviCide1™	Not recommended for use
L25x	CaviWipes1™	Not recommended for use
L25xp	CaviCide1™	Not recommended for use
L25xp	CaviWipes1™	Not recommended for use
L38n	CaviCide1™	Not recommended for use
L38n	CaviWipes1™	Not recommended for use
L38xp	CaviCide1™	Not recommended for use
L52x	CaviCide1™	Not recommended for use
L52x	CaviWipes1™	Not recommended for use
IC10-3	CaviCide1™	Not recommended for use
IC10-3	CaviWipes1™	Not recommended for use
IC10-3	CaviCide1™	Not recommended for use
IC10-3	CaviWipes1™	Not recommended for use
ICTn	CaviCide1™	Not recommended for use
ICTn	CaviWipes1™	Not recommended for use
ICTn	CaviCide1™	Not recommended for use
ICTn	CaviWipes1™	Not recommended for use
P21V	CaviCide1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
P21V	CaviWipes1™	Not recommended for use
iViz System Protective Case	CaviCide™	Not recommended for use
iViz System Protective Case	CaviWipes™	Not recommended for use
iViz System Protective Case	Envirocide™	Not recommended for use
iViz System Protective Case	CaviCide1™	Not recommended for use
iViz System Protective Case	CaviWipes1™	Not recommended for use
iViz System Carry Case	CaviCide1™	Not recommended for use
iViz System Carry Case	CaviWipes1™	Not recommended for use
L12-3	CaviCide1™	Not recommended for use
L12-3	CaviWipes1™	Not recommended for use
L13-6	CaviCide1™	Not recommended for use
L13-6	CaviWipes1™	Not recommended for use
L15-4	CaviCide1™	Not recommended for use
L15-4	CaviWipes1™	Not recommended for use
L19-5	CaviCide1™	Not recommended for use
L19-5	CaviWipes1™	Not recommended for use
L25n	CaviCide1™	Not recommended for use
L38xp	CaviWipes1™	Not recommended for use
L52x Armored	CaviCide™	Not recommended for use
L52x Armored	CaviWipes™	Not recommended for use
L52x Armored	Envirocide™	Not recommended for use
L52x Armored	CaviCide1™	Not recommended for use
L52x Armored	CaviWipes1™	Not recommended for use
L52n	CaviCide1™	Not recommended for use
L52n	CaviWipes1™	Not recommended for use
HSL25x	CaviWipes1™	Not recommended for use



<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
HSL50xp	CaviCide1™	Not recommended for use
HSL50xp	CaviWipes1™	Not recommended for use
HSL25xp	CaviCide1™	Not recommended for use
HSL25xp	CaviWipes1™	Not recommended for use
ICTn	CaviCide1™	Not recommended for use
ICTn	CaviWipes1™	Not recommended for use
FC1 Stand	CaviCide1™	Not recommended for use
FC1 Stand	CaviWipes1™	Not recommended for use
FC1 System	CaviCide™	Not recommended for use
FC1 System	CaviWipes™	Not recommended for use
FC1 System	Envirocide™	Not recommended for use
FC1 System	CaviCide1™	Not recommended for use
FC1 System	CaviWipes1™	Not recommended for use
PowerPack	CaviWipes1™	Not recommended for use
rP19x	CaviCide1™	Not recommended for use
rP19x	CaviWipes1™	Not recommended for use
rP19xp	CaviCide1™	Not recommended for use
rP19xp	CaviWipes1™	Not recommended for use
TTC - S Series	CaviCide1™	Not recommended for use
TTC - S Series	CaviWipes1™	Not recommended for use
SII System	CaviCide1™	Not recommended for use
P11x	CaviWipes1™	Not recommended for use
P21n	CaviCide1™	Not recommended for use
P21n	CaviWipes1™	Not recommended for use
P21x	CaviCide1™	Not recommended for use
P21x	CaviWipes1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
P21xp	CaviCide1™	Not recommended for use
P21xp	CaviWipes1™	Not recommended for use
P5-1	CaviCide1™	Not recommended for use
P5-1	CaviWipes1™	Not recommended for use
PowerPack	CaviCide1™	Not recommended for use
ICTx Revision 80 & below	CaviCide1™	Not recommended for use
ICTx Revision 80 & below	CaviWipes1™	Not recommended for use
ICTx Revision 80 & below	CaviCide1™	Not recommended for use
ICTx Revision 80 & below	CaviWipes1™	Not recommended for use
ICTx Revision 80 & below	CaviCide1™	Not recommended for use
ICTx Revision 80 & below	CaviWipes1™	Not recommended for use
ICTx Revision 81	CaviCide1™	Not recommended for use
ICTx Revision 81	CaviWipes1™	Not recommended for use
ICTx Revision 81	CaviCide1™	Not recommended for use
SII System	CaviWipes1™	Not recommended for use
S-Series TTC	CaviCide1™	Not recommended for use
S-Series TTC	CaviWipes1™	Not recommended for use
V-Universal Stand	CaviCide™	Not recommended for use
V-Universal Stand	CaviWipes™	Not recommended for use
V-Universal Stand	Envirocide™	Not recommended for use
V-Universal Stand	CaviCide1™	Not recommended for use
V-Universal Stand	CaviWipes1™	Not recommended for use
ICTx Revision 81	CaviWipes1™	Not recommended for use
M-Turbo	CaviCide™	Not recommended for use
M-Turbo	CaviWipes™	Not recommended for use
M-Turbo	Envirocide™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
M-Turbo	CaviCide1™	Not recommended for use
M-Turbo	CaviWipes1™	Not recommended for use
NanoMaxx	CaviCide1™	Not recommended for use
NanoMaxx	CaviWipes1™	Not recommended for use
UHF-22	CaviWipes1™	Not recommended for use
UHF-48	CaviCide1™	Not recommended for use
UHF-48	CaviWipes1™	Not recommended for use
UHF-70	CaviCide1™	Not recommended for use
UHF-70	CaviWipes1™	Not recommended for use
P10x	CaviCide1™	Not recommended for use
P10x	CaviWipes1™	Not recommended for use
P10xp	CaviCide1™	Not recommended for use
P10xp	CaviWipes1™	Not recommended for use
P11x	CaviCide1™	Not recommended for use
UHF-22	CaviCide1™	Not recommended for use
EDGE	CaviWipes1™	Not recommended for use
iViz	CaviCide™	Not recommended for use
iViz	CaviWipes™	Not recommended for use
iViz	Envirocide™	Not recommended for use
iViz	CaviCide1™	Not recommended for use
iViz	CaviWipes1™	Not recommended for use
S Series	CaviCide™	Not recommended for use
S Series	CaviWipes™	Not recommended for use
S Series	Envirocide™	Not recommended for use
S Series	CaviCide1™	Not recommended for use
S Series	CaviWipes1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
X-Porte System	CaviCide1™	Not recommended for use
X-Porte System	CaviWipes1™	Not recommended for use
Vevo MD System	CaviCide1™	Not recommended for use
Vevo MD System	CaviWipes1™	Not recommended for use
X-Porte Stand	CaviCide™	Not recommended for use
X-Porte Stand	CaviWipes™	Not recommended for use
X-Porte Stand	Envirocide™	Not recommended for use
X-Porte Stand	CaviCide1™	Not recommended for use
X-Porte Stand	CaviWipes1™	Not recommended for use
EDGE	CaviCide1™	Not recommended for use

<b>GE Healthcare</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
3.0T Signa HDxt	CaviCide™	Not recommended for use
3.0T Signa HDxt	CaviWipes™	Not recommended for use
3.0T Signa HDxt	Envirocide™	Not recommended for use
3.0T Signa HDxt	CaviCide1™	Not recommended for use
3.0T Signa HDxt	CaviWipes1™	Not recommended for use
3.0T Signa HDxt	CaviWipes™ 2.0	Not recommended for use
1.5T Signa HDxt	CaviCide™	Not recommended for use
1.5T Signa HDxt	CaviWipes™	Not recommended for use
1.5T Signa HDxt	Envirocide™	Not recommended for use
1.5T Signa HDxt	CaviCide1™	Not recommended for use
1.5T Signa HDxt	CaviWipes1™	Not recommended for use
1.5T Signa HDxt	CaviWipes™ 2.0	Not recommended for use
Brivo MR355	CaviWipes1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Brivo MR355	CaviWipes™ 2.0	Not recommended for use
Brivo MR355	CaviCide™	Not recommended for use
Brivo MR355	CaviWipes™	Not recommended for use
Brivo MR355	Envirocide™	Not recommended for use
Brivo MR355	CaviCide1™	Not recommended for use
Dash 3000	CaviCide™	Not recommended for use
Dash 3000	CaviWipes™	Not recommended for use
Dash 3000	Envirocide™	Not recommended for use
Dash 3000	CaviCide1™	Not recommended for use
Dash 3000	CaviWipes1™	Not recommended for use
Dash 3000	CaviWipes™ 2.0	Not recommended for use
Dash 4000	CaviCide™	Not recommended for use
Dash 4000	CaviWipes™	Not recommended for use
Dash 4000	Envirocide™	Not recommended for use
Dash 4000	CaviCide1™	Not recommended for use
Dash 4000	CaviWipes1™	Not recommended for use
Dash 4000	CaviWipes™ 2.0	Not recommended for use
Dash 5000	CaviCide™	Not recommended for use
Dash 5000	CaviWipes™	Not recommended for use
Dash 5000	Envirocide™	Not recommended for use
Dash 5000	CaviCide1™	Not recommended for use
Dash 5000	CaviWipes1™	Not recommended for use
Dash 5000	CaviWipes™ 2.0	Not recommended for use
Discovery MR450 1.5T	CaviCide™	Not recommended for use
Discovery MR450 1.5T	CaviWipes™	Not recommended for use
Discovery MR450 1.5T	Envirocide™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Discovery MR450 1.5T	CaviCide1™	Not recommended for use
Discovery MR450 1.5T	CaviWipes1™	Not recommended for use
Discovery MR450 1.5T	CaviWipes™ 2.0	Not recommended for use
Discovery MR750 3.0T	CaviCide™	Not recommended for use
Discovery MR750 3.0T	CaviWipes™	Not recommended for use
Discovery MR750 3.0T	Envirocide™	Not recommended for use
Discovery MR750 3.0T	CaviCide1™	Not recommended for use
Discovery MR750 3.0T	CaviWipes1™	Not recommended for use
Discovery MR750 3.0T	CaviWipes™ 2.0	Not recommended for use
Discovery MR750w 3.0T	CaviCide™	Not recommended for use
Discovery MR750w 3.0T	CaviWipes™	Not recommended for use
Discovery MR750w 3.0T	Envirocide™	Not recommended for use
Discovery MR750w 3.0T	CaviCide1™	Not recommended for use
Discovery MR750w 3.0T	CaviWipes1™	Not recommended for use
Discovery MR750w 3.0T	CaviWipes™ 2.0	Not recommended for use
Flex Coil Positioner	CaviCide™	Not recommended for use
Flex Coil Positioner	CaviWipes™	Not recommended for use
Flex Coil Positioner	Envirocide™	Not recommended for use
Flex Coil Positioner	CaviCide1™	Not recommended for use
Flex Coil Positioner	CaviWipes1™	Not recommended for use
Flex Coil Positioner	CaviWipes™ 2.0	Not recommended for use
LOGIQ S7 R1-R2	CaviCide1™	Not recommended for use
LOGIQ S7 R1-R2	CaviWipes1™	Not recommended for use
LOGIQ S8 R1-R3 (Blue Buttons)	CaviCide1™	Not recommended for use
LOGIQ S8 R1-R3 (Blue Buttons)	CaviWipes1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
LOGIQ S7 R3	CaviCide1™	Not recommended for use
LOGIQ S7 R3	CaviWipes1™	Not recommended for use
LOGIQ S7 XDclr2.0	CaviCide1™	Not recommended for use
LOGIQ S7 XDclr2.0	CaviWipes1™	Not recommended for use
LOGIQ S8R4 Blk Buttons	CaviCide1™	Not recommended for use
LOGIQ S8R4 Blk Buttons	CaviWipes1™	Not recommended for use
LOGIQ S8R4 Gray Buttons	CaviCide1™	Not recommended for use
LOGIQ S8R4 Gray Buttons	CaviWipes1™	Not recommended for use
Magnetic Resonance (MR) Coils	CaviCide™	Not recommended for use
Magnetic Resonance (MR) Coils	CaviWipes™	Not recommended for use
Magnetic Resonance (MR) Coils	Envirocide™	Not recommended for use
Magnetic Resonance (MR) Coils	CaviCide1™	Not recommended for use
Magnetic Resonance (MR) Coils	CaviWipes1™	Not recommended for use
Magnetic Resonance (MR) Coils	CaviWipes™ 2.0	Not recommended for use
MAC 500	CaviWipes™ HP	Not recommended for use
MR Dielectric Pad - Abdomen	CaviCide™	Not recommended for use
MR Dielectric Pad - Abdomen	CaviWipes™	Not recommended for use
MR Dielectric Pad - Abdomen	Envirocide™	Not recommended for use
MR Dielectric Pad - Abdomen	CaviCide1™	Not recommended for use
MR Dielectric Pad - Abdomen	CaviWipes1™	Not recommended for use
MR Dielectric Pad - Abdomen	CaviWipes™ 2.0	Not recommended for use
MR Dielectric Pad - Neck	CaviCide™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
MR Dielectric Pad - Neck	CaviWipes™	Not recommended for use
MR Dielectric Pad - Neck	Envirocide™	Not recommended for use
MR Dielectric Pad - Neck	CaviCide1™	Not recommended for use
MR Dielectric Pad - Neck	CaviWipes1™	Not recommended for use
MR Dielectric Pad - Neck	CaviWipes™ 2.0	Not recommended for use
Optima MR360	CaviWipes1™	Not recommended for use
Optima MR360	CaviWipes™ 2.0	Not recommended for use
Optima MR450w 1.5T	CaviCide™	Not recommended for use
Optima MR450w 1.5T	CaviWipes™	Not recommended for use
Optima MR450w 1.5T	Envirocide™	Not recommended for use
Optima MR450w 1.5T	CaviCide1™	Not recommended for use
Optima MR450w 1.5T	CaviWipes1™	Not recommended for use
Optima MR450w 1.5T	CaviWipes™ 2.0	Not recommended for use
Optima MR360	CaviCide™	Not recommended for use
Optima MR360	CaviWipes™	Not recommended for use
Optima MR360	Envirocide™	Not recommended for use
Optima MR360	CaviCide1™	Not recommended for use
Pediatric Positioner Pad	CaviCide™	Not recommended for use
Pediatric Positioner Pad	CaviWipes™	Not recommended for use
Pediatric Positioner Pad	Envirocide™	Not recommended for use
Pediatric Positioner Pad	CaviCide1™	Not recommended for use
Pediatric Positioner Pad	CaviWipes1™	Not recommended for use
Pediatric Positioner Pad	CaviWipes™ 2.0	Not recommended for use
SIGNA Premier	CaviWipes1™	Not recommended for use
SIGNA Premier	CaviWipes™ 2.0	Not recommended for use
SIGNA Voyager	CaviCide™	Not recommended for use



<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
SIGNA Voyager	CaviWipes™	Not recommended for use
SIGNA Voyager	Envirocide™	Not recommended for use
SIGNA Voyager	CaviCide1™	Not recommended for use
SIGNA Voyager	CaviWipes1™	Not recommended for use
SIGNA Voyager	CaviWipes™ 2.0	Not recommended for use
Tonoport V Blood Pressure system CUFFS	CaviWipes™ HP	Not recommended for use
SIGNA Architect	CaviCide™	Not recommended for use
SIGNA Architect	CaviWipes™	Not recommended for use
SIGNA Architect	Envirocide™	Not recommended for use
SIGNA Architect	CaviCide1™	Not recommended for use
SIGNA Architect	CaviWipes1™	Not recommended for use
SIGNA Architect	CaviWipes™ 2.0	Not recommended for use
SIGNA Artist	CaviCide™	Not recommended for use
SIGNA Artist	CaviWipes™	Not recommended for use
SIGNA Artist	Envirocide™	Not recommended for use
SIGNA Artist	CaviCide1™	Not recommended for use
SIGNA PET / MR	CaviCide™	Not recommended for use
SIGNA PET / MR	CaviWipes™	Not recommended for use
SIGNA PET / MR	Envirocide™	Not recommended for use
SIGNA PET / MR	CaviCide1™	Not recommended for use
SIGNA PET / MR	CaviWipes1™	Not recommended for use
SIGNA PET / MR	CaviWipes™ 2.0	Not recommended for use
SIGNA Pioneer	CaviCide™	Not recommended for use
SIGNA Pioneer	CaviWipes™	Not recommended for use
SIGNA Pioneer	Envirocide™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
SIGNA Pioneer	CaviCide1™	Not recommended for use
SIGNA Pioneer	CaviWipes1™	Not recommended for use
SIGNA Pioneer	CaviWipes™ 2.0	Not recommended for use
SIGNA Premier	CaviCide™	Not recommended for use
SIGNA Premier	CaviWipes™	Not recommended for use
SIGNA Premier	Envirocide™	Not recommended for use
SIGNA Premier	CaviCide1™	Not recommended for use
SIGNA Artist	CaviWipes1™	Not recommended for use
SIGNA Artist	CaviWipes™ 2.0	Not recommended for use
SIGNA Creator	CaviCide™	Not recommended for use
SIGNA Creator	CaviWipes™	Not recommended for use
SIGNA Creator	Envirocide™	Not recommended for use
SIGNA Creator	CaviCide1™	Not recommended for use
SIGNA Creator	CaviWipes1™	Not recommended for use
SIGNA Creator	CaviWipes™ 2.0	Not recommended for use
SIGNA Explorer	CaviCide™	Not recommended for use
SIGNA Explorer	CaviWipes™	Not recommended for use
SIGNA Explorer	Envirocide™	Not recommended for use
SIGNA Explorer	CaviCide1™	Not recommended for use
SIGNA Explorer	CaviWipes1™	Not recommended for use
SIGNA Explorer	CaviWipes™ 2.0	Not recommended for use
Vivid S60	CaviWipes™ Bleach	Not recommended for use
Vivid 70	CaviWipes™ Bleach	Not recommended for use
Vivid 60N	CaviWipes™ Bleach	Not recommended for use
Vivid 70N	CaviWipes™ Bleach	Not recommended for use

### **Hill-Rom**

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Advanta 2 P1190/head & foot boards	CaviCide™	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviWipes™	Not recommended for use
Advanta 2 P1190/head & foot boards	Envirocide™	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviCide1™	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviWipes1™	Not recommended for use
Advanta 2 P1190/head & foot boards	CaviWipes™ 2.0	Not recommended for use

<b>LifeScan</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Onetouch Ping	CaviWipes™ Bleach	Not recommended for use

<b>Marus</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
DC 1700 Series	CaviWipes™ HP	Not recommended for use
DC 1700 Series	CaviWipes™ Bleach	Not recommended for use

<b>Medtronics/ Covidien</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Newport e360S	CaviCide™	Not recommended for use
Newport e360S	CaviWipes™	Not recommended for use
Newport e360S	Envirocide™	Not recommended for use
Newport e360S	CaviCide1™	Not recommended for use

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Newport e360S	CaviWipes1™	Not recommended for use
Newport e360S	CaviWipes™ 2.0	Not recommended for use
Newport e360S	CaviWipes™ Bleach	Not recommended for use
Newport e360P	CaviCide™	Not recommended for use
Newport e360P	CaviWipes™	Not recommended for use
Newport e360P	Envirocide™	Not recommended for use
Newport e360P	CaviCide1™	Not recommended for use
Newport e360P	CaviWipes1™	Not recommended for use
Newport e360P	CaviWipes™ 2.0	Not recommended for use
Newport e360P	CaviWipes™ Bleach	Not recommended for use
Newport e360E	CaviCide™	Not recommended for use
Newport e360E	CaviWipes™	Not recommended for use
Newport e360E	Envirocide™	Not recommended for use
Newport e360E	CaviCide1™	Not recommended for use
Newport e360E	CaviWipes1™	Not recommended for use
Newport e360E	CaviWipes™ 2.0	Not recommended for use
Newport e360E	CaviWipes™ Bleach	Not recommended for use

<b>Midmark</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
True Definition Scanner	CaviCide™	Not recommended for use
True Definition Scanner	CaviWipes™	Not recommended for use
True Definition Scanner	Envirocide™	Not recommended for use

<b>NDS</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>

Radiance 19" Medical Display	CaviCide1™	Not recommended for use
Radiance 19" Medical Display	CaviWipes1™	Not recommended for use
Radiance 19" Medical Display	CaviWipes™ 2.0	Not recommended for use

Philips		
Equipment	Metrex Product	Compatibility Status
Intellivue MP2	CaviWipes™ Bleach	Not recommended for use
MP5T	CaviWipes™ Bleach	Not recommended for use
VS3- Predictive & Tympanic Temperature module-thermometer	CaviWipes™ Bleach	Not recommended for use

Royal		
Equipment	Metrex Product	Compatibility Status
all upholstered chairs & stools	CaviWipes™ Bleach	Not recommended for use
Royal Alliant	CaviWipes™ Bleach	Not recommended for use

SDI		
Equipment	Metrex Product	Compatibility Status
Radii Plus	CaviCide™	Not recommended for use
Radii Plus	Envirocide™	Not recommended for use
Radii Plus	CaviCide1™	Not recommended for use

Siemens		
---------	--	--

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Mammomat Mammography Systems - Compression Trough	CaviWipes™ HP	Not recommended for use
Mammomat Mammography Systems - Standard needle holder, special needle holder for dedicated biopsy systems	CaviWipes™ HP	Not recommended for use
Acuson P500	CaviCide™	Not recommended for use
Acuson P500	CaviWipes™	Not recommended for use
Acuson P500	Envirocide™	Not recommended for use

<b>Sirona</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Schick Sensor	CaviCide1™	Not recommended for use
Schick Sensor	CaviWipes1™	Not recommended for use
XIOS XG sensor & cable	CaviCide1™	Not recommended for use
XIOS XG sensor & cable	CaviWipes1™	Not recommended for use
Orthophos SL 2D	CaviWipes™ HP	Not recommended for use
Orthophos SL 2D	CaviWipes™ Bleach	Not recommended for use
Orthophos SL 3D	CaviWipes™ HP	Not recommended for use
Orthophos SL 3D	CaviWipes™ Bleach	Not recommended for use

<b>Ultradent</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Valo	CaviWipes™ Bleach	Not recommended for use

<b>Welch Allyn</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>

Propaq encore

CaviWipes™ Bleach

Not recommended for use