



Date: December 3, 2012  
To: CaviCide1 Customer  
From: Metrex R&D  
Subject: CaviCide1 and Carcinogenicity

---

CaviCide1 is an Environmental Protection Agency (EPA) registered cleaner and disinfectant. It is used to clean and disinfect inanimate objects, equipment and instrumentation in medical and dental settings. Its active ingredients are isopropyl alcohol, ethanol and quaternary ammonium chloride. A cleaner is contained within the solution.

As stated on the Materials Safety Data Sheet (MSDS) for CaviCide1 Surface Disinfectant/ Decontaminant Cleaner, CaviCide1 is not listed with the Department of Health and Human Services - National Toxicology Program (NTP). CaviCide1 is not listed on the International Agency for Research on Cancer (IARC) Monographs, nor is it listed as an Occupational Safety and Health Administration (OSHA) regulated carcinogen. Additionally, CaviCide1 or any of its components are not listed with the American Conference of Governmental Industrial Hygienists (ACGIH).

Based on this assessment of regulatory agencies governing CaviCide1, neither CaviCide1 nor its components are known carcinogens or potential carcinogens.

If you have additional questions concerning this matter, please contact Metrex Customer Care at 800.841.1428