

Material Compatibility Status

Siemens		
Equipment	Metrex Product	Compatibility Status
ACUSON Juniper	CaviCide™	Compatible
ACUSON Juniper	CaviWipes™	Compatible
ACUSON Juniper	Envirocide™	Compatible
ACUSON Juniper	CaviCide1™	No information
ACUSON Juniper	CaviWipes1™	No information
ACUSON Juniper	CaviWipes™ 2.0	No information
ACUSON Juniper	CaviWipes™ HP	No information
ACUSON Juniper	CaviWipes™ Bleach	No information