

Material Compatibility Status

Siemens		
Equipment	Metrex Product	Compatibility Status
Acuson NX2	CaviCide™	Compatible
Acuson NX2	CaviWipes™	Compatible
Acuson NX2	Envirocide™	Compatible
Acuson NX2	CaviCide1™	Compatible
Acuson NX2	CaviWipes1™	Compatible
Acuson NX2	CaviWipes™ 2.0	Compatible
Acuson NX2	CaviWipes™ HP	No information
Acuson NX2	CaviWipes™ Bleach	Compatible