

Material Compatibility Status

BD		
Equipment	Metrex Product	Compatibility Status
Pyxis Supply Technologies	CaviCide™	Compatible
Pyxis Supply Technologies	CaviWipes™	Compatible
Pyxis Supply Technologies	Envirocide™	Compatible
Pyxis Supply Technologies	CaviCide1™	No information
Pyxis Supply Technologies	CaviWipes1™	No information
Pyxis Supply Technologies	CaviWipes™ 2.0	No information
Pyxis Supply Technologies	CaviWipes™ HP	No information
Pyxis Supply Technologies	CaviWipes™ Bleach	No information
Alaris System	CaviCide™	No information
Alaris System	CaviWipes™	No information
Alaris System	Envirocide™	No information
Alaris System	CaviCide1™	No information
Alaris System	CaviWipes1™	No information
Alaris System	CaviWipes™ 2.0	No information
Alaris System	CaviWipes™ HP	No information
Alaris System	CaviWipes™ Bleach	No information
Pyxis Anesthesia System	CaviCide™	No information
Pyxis Anesthesia System	CaviWipes™	No information
Pyxis Anesthesia System	Envirocide™	No information
Pyxis Anesthesia System	CaviCide1™	No information
Pyxis Anesthesia System	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Pyxis Anesthesia System	CaviWipes™ 2.0	No information
Pyxis Anesthesia System	CaviWipes™ HP	No information
Pyxis Anesthesia System	CaviWipes™ Bleach	No information
Pyxis MedStation	CaviCide™	No information
Pyxis MedStation	CaviWipes™	No information
Pyxis MedStation	Envirocide™	No information
Pyxis MedStation	CaviCide1™	No information
Pyxis MedStation	CaviWipes1™	No information
Pyxis MedStation	CaviWipes™ 2.0	No information
Pyxis MedStation	CaviWipes™ HP	No information
Pyxis MedStation	CaviWipes™ Bleach	No information
Pyxis Supply Technologies	CaviCide™	Compatible
Pyxis Supply Technologies	CaviWipes™	Compatible
Pyxis Supply Technologies	Envirocide™	Compatible
Pyxis Supply Technologies	CaviCide1™	No information
Pyxis Supply Technologies	CaviWipes1™	No information
Pyxis Supply Technologies	CaviWipes™ 2.0	No information
Pyxis Supply Technologies	CaviWipes™ HP	No information
Pyxis Supply Technologies	CaviWipes™ Bleach	No information

Kendall SCD		
Equipment	Metrex Product	Compatibility Status
700 Sequential Compression system	CaviCide™	No information
700 Sequential Compression system	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
700 Sequential Compression system	Envirocide™	No information
700 Sequential Compression system	CaviCide1™	No information
700 Sequential Compression system	CaviWipes1™	No information
700 Sequential Compression system	CaviWipes™ 2.0	No information
700 Sequential Compression system	CaviWipes™ HP	Compatible
700 Sequential Compression system	CaviWipes™ Bleach	Compatible
Model 7325	CaviCide™	No information
Model 7325	CaviWipes™	No information
Model 7325	Envirocide™	No information
Model 7325	CaviCide1™	No information
Model 7325	CaviWipes1™	No information
Model 7325	CaviWipes™ 2.0	No information
Model 7325	CaviWipes™ HP	No information
Model 7325	CaviWipes™ Bleach	No information

SECHRIST Industries		
Equipment	Metrex Product	Compatibility Status
3300 H Monoplace	CaviCide™	No information
3300 H Monoplace	CaviWipes™	No information
3300 H Monoplace	Envirocide™	No information
3300 H Monoplace	CaviCide1™	No information
3300 H Monoplace	CaviWipes1™	No information
3300 H Monoplace	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
3300 H Monoplace	CaviWipes™ HP	Compatible
3300 H Monoplace	CaviWipes™ Bleach	Compatible
4100 HR Monoplace	CaviCide™	No information
4100 HR Monoplace	CaviWipes™	No information
4100 HR Monoplace	Envirocide™	No information
4100 HR Monoplace	CaviCide1™	No information
4100 HR Monoplace	CaviWipes1™	No information
4100 HR Monoplace	CaviWipes™ 2.0	No information
4100 HR Monoplace	CaviWipes™ HP	Compatible
4100 HR Monoplace	CaviWipes™ Bleach	Compatible
Hydraulic Model P/N 21266	CaviCide™	No information
Hydraulic Model P/N 21266	CaviWipes™	No information
Hydraulic Model P/N 21266	Envirocide™	No information
Hydraulic Model P/N 21266	CaviCide1™	No information
Hydraulic Model P/N 21266	CaviWipes1™	No information
Hydraulic Model P/N 21266	CaviWipes™ 2.0	No information
Hydraulic Model P/N 21266	CaviWipes™ HP	No information
Hydraulic Model P/N 21266	CaviWipes™ Bleach	No information
Hydraulic Model P/N 21465	CaviCide™	No information
Hydraulic Model P/N 21465	CaviWipes™	No information
Hydraulic Model P/N 21465	Envirocide™	No information
Hydraulic Model P/N 21465	CaviCide1™	No information
Hydraulic Model P/N 21465	CaviWipes1™	No information
Hydraulic Model P/N 21465	CaviWipes™ 2.0	No information
Hydraulic Model P/N 21465	CaviWipes™ HP	No information
Hydraulic Model P/N 21465	CaviWipes™ Bleach	No information
Standard Model P/N 21464	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Standard Model P/N 21464	CaviWipes™	No information
Standard Model P/N 21464	Envirocide™	No information
Standard Model P/N 21464	CaviCide1™	No information
Standard Model P/N 21464	CaviWipes1™	No information
Standard Model P/N 21464	CaviWipes™ 2.0	No information
Standard Model P/N 21464	CaviWipes™ HP	No information
Standard Model P/N 21464	CaviWipes™ Bleach	No information
Standard Model P/N 21487	CaviCide™	No information
Standard Model P/N 21487	CaviWipes™	No information
Standard Model P/N 21487	Envirocide™	No information
Standard Model P/N 21487	CaviCide1™	No information
Standard Model P/N 21487	CaviWipes1™	No information
Standard Model P/N 21487	CaviWipes™ 2.0	No information
Standard Model P/N 21487	CaviWipes™ HP	No information
Standard Model P/N 21487	CaviWipes™ Bleach	No information
Standard Model P/N 21502	CaviCide™	No information
Standard Model P/N 21502	CaviWipes™	No information
Standard Model P/N 21502	Envirocide™	No information
Standard Model P/N 21502	CaviCide1™	No information
Standard Model P/N 21502	CaviWipes1™	No information
Standard Model P/N 21502	CaviWipes™ 2.0	No information
Standard Model P/N 21502	CaviWipes™ HP	No information
Standard Model P/N 21502	CaviWipes™ Bleach	No information