

Material Compatibility Status

Arizona Healthcare		
Equipment	Metrex Product	Compatibility Status
SpotOn Temp	CaviCide™	No information
SpotOn Temp	CaviWipes™	No information
SpotOn Temp	Envirocide™	No information
SpotOn Temp	CaviCide1™	No information
SpotOn Temp	CaviWipes1™	No information
SpotOn Temp	CaviWipes™ 2.0	No information
SpotOn Temp	CaviWipes™ HP	No information
SpotOn Temp	CaviWipes™ Bleach	No information

Baxter		
Equipment	Metrex Product	Compatibility Status
Spectrum IQ	CaviCide™	No information
Spectrum IQ	CaviWipes™	No information
Spectrum IQ	Envirocide™	No information
Spectrum IQ	CaviCide1™	No information
Spectrum IQ	CaviWipes1™	No information
Spectrum IQ	CaviWipes™ 2.0	No information
Spectrum IQ	CaviWipes™ HP	No information
Spectrum IQ	CaviWipes™ Bleach	Compatible
2M91617A	CaviCide™	No information
2M91617A	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
2M91617A	Envirocide™	No information
2M91617A	CaviCide1™	No information
2M91617A	CaviWipes1™	No information
2M91617A	CaviWipes™ 2.0	No information
2M91617A	CaviWipes™ HP	Compatible
2M91617A	CaviWipes™ Bleach	Compatible
AS 50	CaviCide™	No information
AS 50	CaviWipes™	No information
AS 50	Envirocide™	No information
AS 50	CaviCide1™	No information
Flo-Gard 6200	CaviWipes1™	No information
Flo-Gard 6200	CaviWipes™ 2.0	No information
Flo-Gard 6200	CaviWipes™ HP	No information
Flo-Gard 6200	CaviWipes™ Bleach	Compatible
Flo-Gard 6201	CaviCide™	No information
Flo-Gard 6201	CaviWipes™	No information
Flo-Gard 6201	Envirocide™	No information
Flo-Gard 6201	CaviCide1™	No information
Flo-Gard 6201	CaviWipes1™	No information
Flo-Gard 6201	CaviWipes™ 2.0	No information
Flo-Gard 6201	CaviWipes™ HP	No information
Flo-Gard 6201	CaviWipes™ Bleach	Compatible
Flo-Gard 6301 Dual Channel	CaviCide™	No information
Flo-Gard 6301 Dual Channel	CaviWipes™	No information
Flo-Gard 6301 Dual Channel	Envirocide™	No information
Flo-Gard 6301 Dual Channel	CaviCide1™	No information
Flo-Gard 6301 Dual Channel	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Flo-Gard 6301 Dual Channel	CaviWipes™ 2.0	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ HP	No information
Flo-Gard 6301 Dual Channel	CaviWipes™ Bleach	Compatible
iLED 7	CaviCide™	No information
iLED 7	CaviWipes™	No information
iLED 7	Envirocide™	No information
iLED 7	CaviCide1™	No information
iLED 7	CaviWipes1™	No information
iLED 7	CaviWipes™ 2.0	No information
iLED 7	CaviWipes™ HP	No information
iLED 7	CaviWipes™ Bleach	No information
PCA II Pump	CaviCide™	No information
PCA II Pump	CaviWipes™	No information
PCA II Pump	Envirocide™	No information
PCA II Pump	CaviCide1™	No information
PCA II Pump	CaviWipes1™	No information
PCA II Pump	CaviWipes™ 2.0	No information
PCA II Pump	CaviWipes™ HP	No information
PCA II Pump	CaviWipes™ Bleach	No information
Prismaflex Dialysis Machine	CaviCide™	No information
Prismaflex Dialysis Machine	CaviWipes™	No information
Prismaflex Dialysis Machine	Envirocide™	No information
Prismaflex Dialysis Machine	CaviCide1™	No information
AS 50	CaviWipes1™	No information
AS 50	CaviWipes™ 2.0	No information
AS 50	CaviWipes™ HP	No information
AS 50	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
Auto Syringe AS50	CaviCide™	No information
Auto Syringe AS50	CaviWipes™	No information
Auto Syringe AS50	Envirocide™	No information
Auto Syringe AS50	CaviCide1™	No information
Auto Syringe AS50	CaviWipes1™	No information
Auto Syringe AS50	CaviWipes™ 2.0	No information
Auto Syringe AS50	CaviWipes™ HP	No information
Auto Syringe AS50	CaviWipes™ Bleach	Compatible
Colleague 3	CaviCide™	No information
Colleague 3	CaviWipes™	No information
Colleague 3	Envirocide™	No information
Colleague 3	CaviCide1™	No information
Colleague 3	CaviWipes1™	No information
Colleague 3	CaviWipes™ 2.0	No information
Colleague 3	CaviWipes™ HP	No information
Colleague 3	CaviWipes™ Bleach	Compatible
Colleague CXE	CaviCide™	No information
Colleague CXE	CaviWipes™	No information
Colleague CXE	Envirocide™	No information
Colleague CXE	CaviCide1™	No information
Colleague CXE	CaviWipes1™	No information
Colleague CXE	CaviWipes™ 2.0	No information
Colleague CXE	CaviWipes™ HP	Compatible
Colleague CXE	CaviWipes™ Bleach	Compatible
ExactaMix 2400	CaviCide™	No information
ExactaMix 2400	CaviWipes™	No information
ExactaMix 2400	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
ExactaMix 2400	CaviCide1™	No information
ExactaMix 2400	CaviWipes1™	No information
ExactaMix 2400	CaviWipes™ 2.0	No information
ExactaMix 2400	CaviWipes™ HP	No information
ExactaMix 2400	CaviWipes™ Bleach	No information
Flo-Gard 6200	CaviCide™	No information
Flo-Gard 6200	CaviWipes™	No information
Flo-Gard 6200	Envirocide™	No information
Flo-Gard 6200	CaviCide1™	No information
TruLight 5000	CaviWipes1™	No information
TruLight 5000	CaviWipes™ 2.0	No information
TruLight 5000	CaviWipes™ HP	No information
TruLight 5000	CaviWipes™ Bleach	No information
TherMax	CaviCide™	No information
TherMax	CaviWipes™	No information
TherMax	Envirocide™	No information
TherMax	CaviCide1™	No information
TherMax	CaviWipes1™	No information
TherMax	CaviWipes™ 2.0	No information
TherMax	CaviWipes™ HP	No information
TherMax	CaviWipes™ Bleach	No information
TruLight 5000	CaviCide™	No information
TruLight 5000	CaviWipes™	No information
TruLight 5000	Envirocide™	No information
TruLight 5000	CaviCide1™	No information
Prismaflex Dailysis Machine	CaviWipes1™	No information
Prismaflex Dailysis Machine	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Prismaflex Dailysis Machine	CaviWipes™ HP	No information
Prismaflex Dailysis Machine	CaviWipes™ Bleach	Compatible
PrisMax	CaviCide™	Compatible
PrisMax	CaviWipes™	Compatible
PrisMax	Envirocide™	Compatible
PrisMax	CaviCide1™	No information
PrisMax	CaviWipes1™	No information
PrisMax	CaviWipes™ 2.0	No information
PrisMax	CaviWipes™ HP	Compatible
PrisMax	CaviWipes™ Bleach	Compatible
Sigma spectrum	CaviCide™	Compatible
Sigma spectrum	CaviWipes™	Compatible
Sigma spectrum	Envirocide™	Compatible
Sigma spectrum	CaviCide1™	No information
Sigma spectrum	CaviWipes1™	No information
Sigma spectrum	CaviWipes™ 2.0	No information
Sigma spectrum	CaviWipes™ HP	No information
Sigma spectrum	CaviWipes™ Bleach	Compatible