

Material Compatibility Status

AgaMatrix		
Equipment	Metrex Product	Compatibility Status
Jazz	CaviCide™	Not recommended for use
Jazz	CaviWipes™	Not recommended for use
Jazz	Envirocide™	Not recommended for use
Jazz	CaviCide1™	Not recommended for use
Jazz	CaviWipes1™	Not recommended for use
Jazz	CaviWipes™ 2.0	Not recommended for use
Jazz	CaviWipes™ HP	Not recommended for use
Jazz	CaviWipes™ Bleach	Not recommended for use