

Material Compatibility Status

Air Liquide		
Equipment	Metrex Product	Compatibility Status
Monal T60	CaviCide™	No information
Monal T60	CaviWipes™	No information
Monal T60	Envirocide™	No information
Monal T60	CaviCide1™	No information
Monal T60	CaviWipes1™	No information
Monal T60	CaviWipes™ 2.0	No information
Monal T60	CaviWipes™ HP	No information
Monal T60	CaviWipes™ Bleach	No information

Bayer		
Equipment	Metrex Product	Compatibility Status
CONTOUR USB Meter	CaviWipes1™	No information
CONTOUR USB Meter	CaviWipes™ 2.0	No information
CONTOUR USB Meter	CaviWipes™ HP	No information
CONTOUR USB Meter	CaviWipes™ Bleach	No information
Contour XT	CaviCide™	No information
Contour XT	CaviWipes™	No information
Contour XT	Envirocide™	No information
Contour XT	CaviCide1™	No information
Contour XT	CaviWipes1™	No information
Contour XT	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Contour XT	CaviWipes™ HP	No information
Contour XT	CaviWipes™ Bleach	Compatible
DCA 2000+ 5031C	CaviCide™	No information
DCA 2000+ 5031C	CaviWipes™	No information
DCA 2000+ 5031C	Envirocide™	No information
DCA 2000+ 5031C	CaviCide1™	No information
DCA 2000+ 5031C	CaviWipes1™	No information
DCA 2000+ 5031C	CaviWipes™ 2.0	No information
DCA 2000+ 5031C	CaviWipes™ HP	No information
DCA 2000+ 5031C	CaviWipes™ Bleach	Compatible
Hematek Slide Strainer 4488C	CaviCide™	No information
Hematek Slide Strainer 4488C	CaviWipes™	No information
Hematek Slide Strainer 4488C	Envirocide™	No information
Hematek Slide Strainer 4488C	CaviCide1™	No information
Hematek Slide Strainer 4488C	CaviWipes1™	No information
Hematek Slide Strainer 4488C	CaviWipes™ 2.0	No information
Hematek Slide Strainer 4488C	CaviWipes™ HP	No information
Hematek Slide Strainer 4488C	CaviWipes™ Bleach	No information
Medrad Centargo	CaviCide™	Compatible
Medrad Centargo	CaviWipes™	Compatible
Medrad Centargo	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Medrad Centargo	CaviCide1™	Compatible
Medrad Centargo	CaviWipes1™	Compatible
Medrad Centargo	CaviWipes™ 2.0	Compatible
Medrad Centargo	CaviWipes™ HP	No information
Medrad Centargo	CaviWipes™ Bleach	No information
Medrad Mark 7 Arterion	CaviCide™	No information
Medrad Mark 7 Arterion	CaviWipes™	No information
Medrad Mark 7 Arterion	Envirocide™	No information
Medrad Mark 7 Arterion	CaviCide1™	No information
Breeze 2	CaviCide™	No information
Breeze 2	CaviWipes™	No information
Breeze 2	Envirocide™	No information
Breeze 2	CaviCide1™	No information
Breeze 2	CaviWipes1™	No information
Breeze 2	CaviWipes™ 2.0	No information
Breeze 2	CaviWipes™ HP	No information
Breeze 2	CaviWipes™ Bleach	Compatible
Contour	CaviCide™	No information
Contour	CaviWipes™	No information
Contour	Envirocide™	No information
Contour	CaviCide1™	No information
Contour	CaviWipes1™	No information
Contour	CaviWipes™ 2.0	No information
Contour	CaviWipes™ HP	No information
Contour	CaviWipes™ Bleach	Compatible
Contour Link	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Contour Link	CaviWipes™	No information
Contour Link	Envirocide™	No information
Contour Link	CaviCide1™	No information
Contour Link	CaviWipes1™	No information
Contour Link	CaviWipes™ 2.0	No information
Contour Link	CaviWipes™ HP	No information
Contour Link	CaviWipes™ Bleach	No information
Contour Next EZ	CaviCide™	No information
Contour Next EZ	CaviWipes™	No information
Contour Next EZ	Envirocide™	No information
Contour Next EZ	CaviCide1™	No information
Contour Next EZ	CaviWipes1™	No information
Contour Next EZ	CaviWipes™ 2.0	No information
Contour Next EZ	CaviWipes™ HP	No information
Contour Next EZ	CaviWipes™ Bleach	Compatible
CONTOUR USB Meter	CaviCide™	No information
CONTOUR USB Meter	CaviWipes™	No information
CONTOUR USB Meter	Envirocide™	No information
CONTOUR USB Meter	CaviCide1™	No information
Medrad Spectris Solaris EP	CaviCide™	Compatible
Medrad Spectris Solaris EP	CaviWipes™	Compatible
Medrad Spectris Solaris EP	Envirocide™	Compatible
Medrad Spectris Solaris EP	CaviCide1™	No information
Medrad Spectris Solaris EP	CaviWipes1™	No information
Medrad Spectris Solaris EP	CaviWipes™ 2.0	Compatible
Medrad Spectris Solaris EP	CaviWipes™ HP	No information

Equipment	Metrex Product	Compatibility Status
Medrad Spectris Solaris EP	CaviWipes™ Bleach	No information
Stellant ISI 800 System	CaviCide™	No information
Stellant ISI 800 System	CaviWipes™	No information
Stellant ISI 800 System	Envirocide™	No information
Stellant ISI 800 System	CaviCide1™	No information
Stellant ISI 800 System	CaviWipes1™	No information
Stellant ISI 800 System	CaviWipes™ 2.0	No information
Stellant ISI 800 System	CaviWipes™ HP	No information
Stellant ISI 800 System	CaviWipes™ Bleach	No information
Medrad Mark 7 Arterion	CaviWipes1™	No information
Medrad Mark 7 Arterion	CaviWipes™ 2.0	No information
Medrad Mark 7 Arterion	CaviWipes™ HP	No information
Medrad Mark 7 Arterion	CaviWipes™ Bleach	Compatible
Medrad MRX Experian	CaviCide™	No information
Medrad MRX Experian	CaviWipes™	No information
Medrad MRX Experian	Envirocide™	No information
Medrad MRX Experian	CaviCide1™	No information
Medrad MRX Experian	CaviWipes1™	No information
Medrad MRX Experian	CaviWipes™ 2.0	No information
Medrad MRX Experian	CaviWipes™ HP	No information
Medrad MRX Experian	CaviWipes™ Bleach	No information
Medrad MRXperion	CaviCide™	Compatible
Medrad MRXperion	CaviWipes™	Compatible
Medrad MRXperion	Envirocide™	Compatible
Medrad MRXperion	CaviCide1™	No information
Medrad MRXperion	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Medrad MRXperion	CaviWipes™ 2.0	No information
Medrad MRXperion	CaviWipes™ HP	No information
Medrad MRXperion	CaviWipes™ Bleach	No information

Spacelabs		
Equipment	Metrex Product	Compatibility Status
Tri-Cuff	CaviCide™	No information
Tri-Cuff	CaviWipes™	No information
Tri-Cuff	Envirocide™	No information
Tri-Cuff	CaviCide1™	No information
Tri-Cuff	CaviWipes1™	No information
Tri-Cuff	CaviWipes™ 2.0	No information
Tri-Cuff	CaviWipes™ HP	No information
Tri-Cuff	CaviWipes™ Bleach	No information
TRU CUFF	CaviCide™	No information
TRU CUFF	CaviWipes™	No information
TRU CUFF	Envirocide™	No information
TRU CUFF	CaviCide1™	No information
TRU CUFF	CaviWipes1™	No information
TRU CUFF	CaviWipes™ 2.0	No information
TRU CUFF	CaviWipes™ HP	No information
TRU CUFF	CaviWipes™ Bleach	Compatible
TRU LINK	CaviCide™	No information
TRU LINK	CaviWipes™	No information
TRU LINK	Envirocide™	No information
TRU LINK	CaviCide1™	No information
TRU LINK	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
TRU LINK	CaviWipes™ 2.0	No information
TRU LINK	CaviWipes™ HP	No information
TRU LINK	CaviWipes™ Bleach	Compatible
Ultaview SL	CaviCide™	No information
Ultaview SL	CaviWipes™	No information
Ultaview SL	Envirocide™	No information
Ultaview SL	CaviCide1™	No information
Ultaview SL	CaviWipes1™	No information
Ultaview SL	CaviWipes™ 2.0	No information
Ultaview SL	CaviWipes™ HP	No information
Ultaview SL	CaviWipes™ Bleach	Compatible
Ultraview Care Network/Transducers	CaviWipes1™	No information
Ultraview Care Network/Transducers	CaviWipes™ 2.0	No information
Ultraview Care Network/Transducers	CaviWipes™ HP	No information
Ultraview Care Network/Transducers	CaviWipes™ Bleach	No information
Ultraview Care Network/Sensors	CaviCide™	No information
Ultraview Care Network/Sensors	CaviWipes™	No information
Ultraview Care Network/Sensors	Envirocide™	No information
Ultraview Care Network/Sensors	CaviCide1™	No information
Ultraview Care Network/Sensors	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Ultraview Care Network/Sensors	CaviWipes™ 2.0	No information
Ultraview Care Network/Sensors	CaviWipes™ HP	No information
Ultraview Care Network/Sensors	CaviWipes™ Bleach	Compatible
Ultraview Care Network	CaviCide™	No information
Ultraview Care Network	CaviWipes™	No information
Ultraview Care Network	Envirocide™	No information
Ultraview Care Network	CaviCide1™	No information
Ultraview Care Network	CaviWipes1™	No information
Ultraview Care Network	CaviWipes™ 2.0	No information
Ultraview Care Network	CaviWipes™ HP	No information
Ultraview Care Network	CaviWipes™ Bleach	No information
Ultraview Care Network/Monitor	CaviCide™	No information
Ultraview Care Network/Monitor	CaviWipes™	No information
Ultraview Care Network/Monitor	Envirocide™	No information
Ultraview Care Network/Monitor	CaviCide1™	No information
Ultraview Care Network/Monitor	CaviWipes1™	No information
Ultraview Care Network/Monitor	CaviWipes™ 2.0	No information
Ultraview Care Network/Monitor	CaviWipes™ HP	No information
Ultraview Care Network/Monitor	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Ultraview Care Network/Transducers	CaviCide™	No information
Ultraview Care Network/Transducers	CaviWipes™	No information
Ultraview Care Network/Transducers	Envirocide™	No information
Ultraview Care Network/Transducers	CaviCide1™	No information