

# Material Compatibility Status

| ArjoHuntLeigh                             |                |                      |
|-------------------------------------------|----------------|----------------------|
| Equipment                                 | Metrex Product | Compatibility Status |
| Sara 3000                                 | CaviWipes™ 2.0 | No information       |
| Sara-Plus                                 | CaviWipes™ 2.0 | No information       |
| Wipeable Standing sling                   | CaviWipes™ 2.0 | Compatible           |
| Alenti                                    | CaviWipes™ 2.0 | No information       |
| Bari Breeze                               | CaviWipes™ 2.0 | No information       |
| Bari Breeze                               | CaviWipes™ 2.0 | No information       |
| Flowtron Universal                        | CaviWipes™ 2.0 | No information       |
| Intermittent pneumatic compression system | CaviWipes™ 2.0 | No information       |
| Maxi Air                                  | CaviWipes™ 2.0 | No information       |
| Maxi Lite                                 | CaviWipes™ 2.0 | No information       |
| Maxi Sky 1000                             | CaviWipes™ 2.0 | No information       |
| Flowtron ACS800 + Tri Pulse               | CaviWipes™ 2.0 | No information       |
| Flowtron ACS900                           | CaviWipes™ 2.0 | No information       |
| Flowtron Excel                            | CaviWipes™ 2.0 | No information       |
| Maxi sky 1000 Charging Station            | CaviWipes™ 2.0 | No information       |
| Maxi Twin Lift                            | CaviWipes™ 2.0 | No information       |

| graham-field/Hausted |                |                      |
|----------------------|----------------|----------------------|
| Equipment            | Metrex Product | Compatibility Status |

|                                                     |                |                |
|-----------------------------------------------------|----------------|----------------|
| ESC2 Powered<br>Ophthalmology Platform              | CaviWipes™ 2.0 | Compatible     |
| Fluoro-Track Fluoroscopy<br>Stretcher Model: 493RPC | CaviWipes™ 2.0 | Compatible     |
| Horizon™ AirGlide™<br>Stretcher                     | CaviWipes™ 2.0 | Compatible     |
| Horizon™ Retracto™ Rail<br>Stretcher                | CaviWipes™ 2.0 | Compatible     |
| HSS500 Surgical Stool                               | CaviWipes™ 2.0 | Compatible     |
| 5E82 Electric Surgi-Stretcher<br>Series             | CaviWipes™ 2.0 | Compatible     |
| EPC500ST Hausted<br>PROCEDURE CHAIR                 | CaviWipes™ 2.0 | No information |
| MBC MAMMOGRAPHY /<br>BIOPSY CHAIR                   | CaviWipes™ 2.0 | Compatible     |
| Mobile Hydraulic Surgi-<br>Stretcher (578)          | CaviWipes™ 2.0 | Compatible     |
| Mobile Powered Surgi-<br>Stretcher (5E8)            | CaviWipes™ 2.0 | Compatible     |
| Pediatric Transport Stretcher<br>Model 4160         | CaviWipes™ 2.0 | Compatible     |
| VIC429 VIDEO IMAGING<br>CHAIR                       | CaviWipes™ 2.0 | Compatible     |
| Youth Transport Stretcher                           | CaviWipes™ 2.0 | Compatible     |