

Material Compatibility Status

Armeda		
Equipment	Metrex Product	Compatibility Status
Platinum Breast Multi-user Hospital Grade Electric Breast Pump	CaviCide™	Compatible
Platinum Breast Multi-user Hospital Grade Electric Breast Pump	CaviWipes™	Compatible
Platinum Breast Multi-user Hospital Grade Electric Breast Pump	Envirocide™	Compatible
Platinum Breast Multi-user Hospital Grade Electric Breast Pump	CaviCide1™	No information
Platinum Breast Multi-user Hospital Grade Electric Breast Pump	CaviWipes1™	No information
Platinum Breast Multi-user Hospital Grade Electric Breast Pump	CaviWipes™ 2.0	No information
Platinum Breast Multi-user Hospital Grade Electric Breast Pump	CaviWipes™ HP	No information
Platinum Breast Multi-user Hospital Grade Electric Breast Pump	CaviWipes™ Bleach	No information