

Material Compatibility Status

Cardiac science		
Equipment	Metrex Product	Compatibility Status
Power Heart AED G3 9300E	CaviCide™	No information
Power Heart AED G3 9300E	CaviWipes™	No information
Power Heart AED G3 9300E	Envirocide™	No information
Power Heart AED G3 9300E	CaviCide1™	No information
Power Heart AED G3 9300E	CaviWipes1™	No information
Power Heart AED G3 9300E	CaviWipes™ 2.0	No information
Power Heart AED G3 9300E	CaviWipes™ HP	No information
Power Heart AED G3 9300E	CaviWipes™ Bleach	No information
Quinton Eclipse Premier ECG	CaviCide™	No information
Quinton Eclipse Premier ECG	CaviWipes™	No information
Quinton Eclipse Premier ECG	Envirocide™	No information
Quinton Eclipse Premier ECG	CaviCide1™	No information
Quinton Eclipse Premier ECG	CaviWipes1™	No information
Quinton Eclipse Premier ECG	CaviWipes™ 2.0	No information
Quinton Eclipse Premier ECG	CaviWipes™ HP	No information
Quinton Eclipse Premier ECG	CaviWipes™ Bleach	No information

Fresenius		
Equipment	Metrex Product	Compatibility Status
CRM III Outer Casing	CaviWipes1™	No information
CRM III Outer Casing	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
CRM III Outer Casing	CaviWipes™ HP	No information
CRM III Outer Casing	CaviWipes™ Bleach	No information
GranuFlo Concentrate Dissolution Unit II	CaviCide™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™	No information
GranuFlo Concentrate Dissolution Unit II	Envirocide™	No information
GranuFlo Concentrate Dissolution Unit II	CaviCide1™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes1™	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ 2.0	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ HP	No information
GranuFlo Concentrate Dissolution Unit II	CaviWipes™ Bleach	Compatible
Liberty Cyclor	CaviCide™	Not recommended for use
Liberty Cyclor	CaviWipes™	Not recommended for use
Liberty Cyclor	Envirocide™	Not recommended for use
Liberty Cyclor	CaviCide1™	Not recommended for use
Liberty Cyclor	CaviWipes1™	Not recommended for use
Liberty Cyclor	CaviWipes™ 2.0	Not recommended for use
Liberty Cyclor	CaviWipes™ HP	No information
Liberty Cyclor	CaviWipes™ Bleach	Not recommended for use
Newton Cyclor	CaviCide™	Not recommended for use
Newton Cyclor	CaviWipes™	Not recommended for use
Newton Cyclor	Envirocide™	Not recommended for use
Newton Cyclor	CaviCide1™	Not recommended for use

Equipment	Metrex Product	Compatibility Status
Newton Cyclor	CaviWipes1™	Not recommended for use
Newton Cyclor	CaviWipes™ 2.0	Not recommended for use
Newton Cyclor	CaviWipes™ HP	No information
Newton Cyclor	CaviWipes™ Bleach	Not recommended for use
2008K	CaviCide™	Not recommended for use
2008K	CaviWipes™	Not recommended for use
2008K	Envirocide™	Not recommended for use
2008K	CaviCide1™	Not recommended for use
2008K	CaviWipes1™	Not recommended for use
2008K	CaviWipes™ 2.0	Not recommended for use
2008K	CaviWipes™ HP	No information
2008K	CaviWipes™ Bleach	Not recommended for use
Crit-Line Monitor III Display	CaviCide™	No information
Crit-Line Monitor III Display	CaviWipes™	No information
Crit-Line Monitor III Display	Envirocide™	No information
Crit-Line Monitor III Display	CaviCide1™	No information
Crit-Line Monitor III Display	CaviWipes1™	No information
Crit-Line Monitor III Display	CaviWipes™ 2.0	No information
Crit-Line Monitor III Display	CaviWipes™ HP	No information
Crit-Line Monitor III Display	CaviWipes™ Bleach	Not recommended for use
Crit-Line Monitor III Sensor clip	CaviCide™	Compatible
Crit-Line Monitor III Sensor clip	CaviWipes™	Compatible
Crit-Line Monitor III Sensor clip	Envirocide™	Compatible
Crit-Line Monitor III Sensor clip	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Crit-Line Monitor III Sensor clip	CaviWipes1™	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ 2.0	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ HP	No information
Crit-Line Monitor III Sensor clip	CaviWipes™ Bleach	Compatible
CRM III Outer Casing	CaviCide™	Compatible
CRM III Outer Casing	CaviWipes™	Compatible
CRM III Outer Casing	Envirocide™	Compatible
CRM III Outer Casing	CaviCide1™	No information