

Material Compatibility Status

Covidien		
Equipment	Metrex Product	Compatibility Status
Achieva Portable Ventilator	CaviWipes1™	No information
Achieva Portable Ventilator	CaviWipes™ 2.0	No information
Achieva Portable Ventilator	CaviWipes™ HP	No information
Achieva Portable Ventilator	CaviWipes™ Bleach	No information
BISX/ Display	CaviCide™	No information
BISX/ Display	CaviWipes™	No information
BISX/ Display	Envirocide™	No information
BISX/ Display	CaviCide1™	No information
BISX/ Display	CaviWipes1™	No information
BISX/ Display	CaviWipes™ 2.0	No information
BISX/ Display	CaviWipes™ HP	No information
BISX/ Display	CaviWipes™ Bleach	Compatible
CapnoStream 20P	CaviCide™	No information
CapnoStream 20P	CaviWipes™	No information
CapnoStream 20P	Envirocide™	No information
CapnoStream 20P	CaviCide1™	No information
CapnoStream 20P	CaviWipes1™	No information
CapnoStream 20P	CaviWipes™ 2.0	No information
CapnoStream 20P	CaviWipes™ HP	No information
CapnoStream 20P	CaviWipes™ Bleach	Compatible
Infant PMR2 Recusitator	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
Infant PMR2 Recusitator	CaviWipes™	No information
Infant PMR2 Recusitator	Envirocide™	No information
Infant PMR2 Recusitator	CaviCide1™	No information
Infant PMR2 Recusitator	CaviWipes1™	No information
Infant PMR2 Recusitator	CaviWipes™ 2.0	No information
Infant PMR2 Recusitator	CaviWipes™ HP	No information
Infant PMR2 Recusitator	CaviWipes™ Bleach	No information
INVOS Model 5100C	CaviCide™	No information
INVOS Model 5100C	CaviWipes™	No information
INVOS Model 5100C	Envirocide™	No information
INVOS Model 5100C	CaviCide1™	No information
INVOS Model 5100C	CaviWipes1™	No information
INVOS Model 5100C	CaviWipes™ 2.0	No information
INVOS Model 5100C	CaviWipes™ HP	No information
INVOS Model 5100C	CaviWipes™ Bleach	No information
Kangaroo ePump	CaviCide™	No information
Kangaroo ePump	CaviWipes™	No information
Kangaroo ePump	Envirocide™	No information
Kangaroo ePump	CaviCide1™	No information
Model 560 Portable Ventilator	CaviCide™	No information
Model 560 Portable Ventilator	CaviWipes™	No information
Model 560 Portable Ventilator	Envirocide™	No information
Model 560 Portable Ventilator	CaviCide1™	No information

Equipment	Metrex Product	Compatibility Status
Model 560 Portable Ventilator	CaviWipes1™	No information
Model 560 Portable Ventilator	CaviWipes™ 2.0	No information
Model 560 Portable Ventilator	CaviWipes™ HP	Compatible
Model 560 Portable Ventilator	CaviWipes™ Bleach	Compatible
Nellcor OxiMax N560	CaviCide™	No information
Nellcor OxiMax N560	CaviWipes™	No information
Nellcor OxiMax N560	Envirocide™	No information
Nellcor OxiMax N560	CaviCide1™	No information
Nellcor OxiMax N560	CaviWipes1™	No information
Nellcor OxiMax N560	CaviWipes™ 2.0	No information
Nellcor OxiMax N560	CaviWipes™ HP	No information
Nellcor OxiMax N560	CaviWipes™ Bleach	Compatible
Nellcor PM1000N	CaviCide™	No information
Nellcor PM1000N	CaviWipes™	No information
Nellcor PM1000N	Envirocide™	No information
Nellcor PM1000N	CaviCide1™	No information
Nellcor PM1000N	CaviWipes1™	No information
Nellcor PM1000N	CaviWipes™ 2.0	No information
Nellcor PM1000N	CaviWipes™ HP	No information
Nellcor PM1000N	CaviWipes™ Bleach	Compatible
Nellcor PM1000N	CaviCide™	No information
Nellcor PM1000N	CaviWipes™	No information
Nellcor PM1000N	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Nellcor PM1000N	CaviCide1™	No information
Nellcor PM1000N	CaviWipes1™	No information
Nellcor PM1000N	CaviWipes™ 2.0	No information
Nellcor PM1000N	CaviWipes™ HP	No information
Nellcor PM1000N	CaviWipes™ Bleach	No information
Newport e360/exhalation flow sensor	CaviCide™	No information
Newport e360/exhalation flow sensor	CaviWipes™	No information
Newport e360/exhalation flow sensor	Envirocide™	No information
Newport e360/exhalation flow sensor	CaviCide1™	No information
Newport e360/exhalation flow sensor	CaviWipes1™	No information
Newport e360/exhalation flow sensor	CaviWipes™ 2.0	No information
Newport e360/exhalation flow sensor	CaviWipes™ HP	No information
Newport e360/exhalation flow sensor	CaviWipes™ Bleach	No information
Newport e360 exterior	CaviCide™	No information
Newport e360 exterior	CaviWipes™	No information
Newport e360 exterior	Envirocide™	No information
Newport e360 exterior	CaviCide1™	No information
Newport e360 exterior	CaviWipes1™	No information
Newport e360 exterior	CaviWipes™ 2.0	No information
Newport e360 exterior	CaviWipes™ HP	No information
Newport e360 exterior	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Newport e360/fan filter	CaviCide™	No information
Newport e360/fan filter	CaviWipes™	No information
Newport e360/fan filter	Envirocide™	No information
Newport e360/fan filter	CaviCide1™	No information
Newport e360/fan filter	CaviWipes1™	No information
Newport e360/fan filter	CaviWipes™ 2.0	No information
Newport e360/fan filter	CaviWipes™ HP	No information
Newport e360/fan filter	CaviWipes™ Bleach	No information
Newport e360/sensor cable	CaviCide™	No information
Newport e360/sensor cable	CaviWipes™	No information
Newport e360/sensor cable	Envirocide™	No information
Newport e360/sensor cable	CaviCide1™	No information
Newport e360/sensor cable	CaviWipes1™	No information
Newport e360/sensor cable	CaviWipes™ 2.0	No information
Newport e360/sensor cable	CaviWipes™ HP	No information
Newport e360/sensor cable	CaviWipes™ Bleach	No information
Newport HT70	CaviCide™	No information
Newport HT70	CaviWipes™	No information
Newport HT70	Envirocide™	No information
Newport HT70	CaviCide1™	No information
Newport HT70	CaviWipes1™	No information
Newport HT70	CaviWipes™ 2.0	No information
Newport HT70	CaviWipes™ HP	No information
Newport HT70	CaviWipes™ Bleach	No information
OxiMax N-600x	CaviCide™	No information
OxiMax N-600x	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
OxiMax N-600x	Envirocide™	No information
OxiMax N-600x	CaviCide1™	No information
OxiMax N-600x	CaviWipes1™	No information
OxiMax N-600x	CaviWipes™ 2.0	No information
OxiMax N-600x	CaviWipes™ HP	No information
OxiMax N-600x	CaviWipes™ Bleach	Compatible
Achieva Portable Ventilator	CaviCide™	No information
Achieva Portable Ventilator	CaviWipes™	No information
Achieva Portable Ventilator	Envirocide™	No information
Achieva Portable Ventilator	CaviCide1™	No information
SCD Express	CaviCide™	No information
SCD Express	CaviWipes™	No information
SCD Express	Envirocide™	No information
SCD Express	CaviCide1™	No information
SCD Express	CaviWipes1™	No information
SCD Express	CaviWipes™ 2.0	No information
SCD Express	CaviWipes™ HP	No information
SCD Express	CaviWipes™ Bleach	No information
AV Impulse 6060 Controller & Accessories	CaviCide™	No information
AV Impulse 6060 Controller & Accessories	CaviWipes™	No information
AV Impulse 6060 Controller & Accessories	Envirocide™	No information
AV Impulse 6060 Controller & Accessories	CaviCide1™	No information
AV Impulse 6060 Controller & Accessories	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
AV Impulse 6060 Controller & Accessories	CaviWipes™ 2.0	No information
AV Impulse 6060 Controller & Accessories	CaviWipes™ HP	No information
AV Impulse 6060 Controller & Accessories	CaviWipes™ Bleach	No information
Puritan-Bennett 540 series	CaviCide™	Compatible
Puritan-Bennett 540 series	CaviWipes™	Compatible
Puritan-Bennett 540 series	Envirocide™	Compatible
Puritan-Bennett 540 series	CaviCide1™	No information
Kangaroo ePump	CaviWipes1™	No information
Kangaroo ePump	CaviWipes™ 2.0	No information
Kangaroo ePump	CaviWipes™ HP	No information
Kangaroo ePump	CaviWipes™ Bleach	Compatible
Kendall SCD 700 System	CaviCide™	No information
Kendall SCD 700 System	CaviWipes™	No information
Kendall SCD 700 System	Envirocide™	No information
Kendall SCD 700 System	CaviCide1™	No information
Kendall SCD 700 System	CaviWipes1™	No information
Kendall SCD 700 System	CaviWipes™ 2.0	No information
Kendall SCD 700 System	CaviWipes™ HP	No information
Kendall SCD 700 System	CaviWipes™ Bleach	Compatible
Model 520 Portable Ventilator	CaviCide™	No information
Model 520 Portable Ventilator	CaviWipes™	No information
Model 520 Portable Ventilator	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Model 520 Portable Ventilator	CaviCide1™	No information
Model 520 Portable Ventilator	CaviWipes1™	No information
Model 520 Portable Ventilator	CaviWipes™ 2.0	No information
Model 520 Portable Ventilator	CaviWipes™ HP	Compatible
Model 520 Portable Ventilator	CaviWipes™ Bleach	Compatible
Puritan-Bennet 980	CaviWipes1™	No information
Puritan-Bennet 980	CaviWipes™ 2.0	No information
Puritan-Bennet 980	CaviWipes™ HP	Compatible
Puritan-Bennet 980	CaviWipes™ Bleach	Compatible
Puritan-Bennett 700 Series	CaviCide™	No information
Puritan-Bennett 700 Series	CaviWipes™	No information
Puritan-Bennett 700 Series	Envirocide™	No information
Puritan-Bennett 700 Series	CaviCide1™	No information
Puritan-Bennett 700 Series	CaviWipes1™	No information
Puritan-Bennett 700 Series	CaviWipes™ 2.0	No information
Puritan-Bennett 700 Series	CaviWipes™ HP	No information
Puritan-Bennett 700 Series	CaviWipes™ Bleach	No information
Puritan-Bennett 800 Series/exterior	CaviCide™	Compatible
Puritan-Bennett 800 Series/exterior	CaviWipes™	Compatible
Puritan-Bennett 800 Series/exterior	Envirocide™	Compatible

Equipment	Metrex Product	Compatibility Status
Puritan-Bennett 800 Series/exterior	CaviCide1™	No information
Puritan-Bennett 800 Series/exterior	CaviWipes1™	No information
Puritan-Bennett 800 Series/exterior	CaviWipes™ 2.0	No information
Puritan-Bennett 800 Series/exterior	CaviWipes™ HP	Compatible
Puritan-Bennett 800 Series/exterior	CaviWipes™ Bleach	Compatible
OxiMax N-65	CaviCide™	No information
OxiMax N-65	CaviWipes™	No information
OxiMax N-65	Envirocide™	No information
OxiMax N-65	CaviCide1™	No information
OxiMax N-65	CaviWipes1™	No information
OxiMax N-65	CaviWipes™ 2.0	No information
OxiMax N-65	CaviWipes™ HP	No information
OxiMax N-65	CaviWipes™ Bleach	Compatible
Puritan-Bennet 840	CaviCide™	Compatible
Puritan-Bennet 840	CaviWipes™	Compatible
Puritan-Bennet 840	Envirocide™	Compatible
Puritan-Bennet 840	CaviCide1™	No information
Puritan-Bennet 840	CaviWipes1™	No information
Puritan-Bennet 840	CaviWipes™ 2.0	No information
Puritan-Bennet 840	CaviWipes™ HP	Compatible
Puritan-Bennet 840	CaviWipes™ Bleach	Compatible
Puritan-Bennet 980	CaviCide™	Compatible
Puritan-Bennet 980	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
Puritan-Bennet 980	Envirocide™	Compatible
Puritan-Bennet 980	CaviCide1™	No information
Puritan-Bennett 540 series	CaviWipes1™	No information
Puritan-Bennett 540 series	CaviWipes™ 2.0	No information
Puritan-Bennett 540 series	CaviWipes™ HP	Compatible
Puritan-Bennett 540 series	CaviWipes™ Bleach	Compatible

Physio-Control		
Equipment	Metrex Product	Compatibility Status
LifePak 12	CaviCide™	Compatible
LifePak 12	CaviWipes™	Compatible
LifePak 12	Envirocide™	Compatible
LifePak 12	CaviCide1™	Compatible
LifePak 12	CaviWipes1™	Compatible
LifePak 12	CaviWipes™ 2.0	Compatible
LifePak 12	CaviWipes™ HP	No information
LifePak 12	CaviWipes™ Bleach	No information
LifePak 15	CaviCide™	Compatible
LifePak 15	CaviWipes™	Compatible
LifePak 15	Envirocide™	Compatible
LifePak 15	CaviCide1™	Compatible
LifePak 15	CaviWipes1™	Compatible
LifePak 15	CaviWipes™ 2.0	Compatible
LifePak 15	CaviWipes™ HP	No information
LifePak 15	CaviWipes™ Bleach	No information
LifePak CR Plus Case Only Disposable Pads	CaviCide™	Compatible

Equipment	Metrex Product	Compatibility Status
LifePak CR Plus Case Only Disposable Pads	CaviWipes™	Compatible
LifePak CR Plus Case Only Disposable Pads	Envirocide™	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviCide1™	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviWipes1™	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviWipes™ 2.0	Compatible
LifePak CR Plus Case Only Disposable Pads	CaviWipes™ HP	No information
LifePak CR Plus Case Only Disposable Pads	CaviWipes™ Bleach	No information
LifePak Express Case only disposable pads	CaviCide™	Compatible
LifePak Express Case only disposable pads	CaviWipes™	Compatible
LifePak Express Case only disposable pads	Envirocide™	Compatible
LifePak Express Case only disposable pads	CaviCide1™	Compatible
LifePak Express Case only disposable pads	CaviWipes1™	Compatible
LifePak Express Case only disposable pads	CaviWipes™ 2.0	Compatible
LifePak Express Case only disposable pads	CaviWipes™ HP	No information
LifePak Express Case only disposable pads	CaviWipes™ Bleach	No information