

Material Compatibility Status

Hologic		
Equipment	Metrex Product	Compatibility Status
Atec Sapphire	CaviCide™	No information
Atec Sapphire	CaviWipes™	No information
Atec Sapphire	Envirocide™	No information
Atec Sapphire	CaviCide1™	No information
Atec Sapphire	CaviWipes1™	No information
Atec Sapphire	CaviWipes™ 2.0	No information
Atec Sapphire	CaviWipes™ HP	No information
Atec Sapphire	CaviWipes™ Bleach	No information
Brevera	CaviCide™	No information
Brevera	CaviWipes™	No information
Brevera	Envirocide™	No information
Brevera	CaviCide1™	No information
Brevera	CaviWipes1™	No information
Brevera	CaviWipes™ 2.0	No information
Brevera	CaviWipes™ HP	No information
Brevera	CaviWipes™ Bleach	Compatible
Faxitron CT	CaviCide™	No information
Faxitron CT	CaviWipes™	No information
Faxitron CT	Envirocide™	No information
Faxitron CT	CaviCide1™	No information
Faxitron CT	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Faxitron CT	CaviWipes™ 2.0	No information
Faxitron CT	CaviWipes™ HP	Compatible
Faxitron CT	CaviWipes™ Bleach	Compatible
Faxitron Path Table	CaviCide™	No information
Faxitron Path Table	CaviWipes™	No information
Faxitron Path Table	Envirocide™	No information
Faxitron Path Table	CaviCide1™	No information
Faxitron Path Table	CaviWipes1™	No information
Faxitron Path Table	CaviWipes™ 2.0	No information
Faxitron Path Table	CaviWipes™ HP	No information
Faxitron Path Table	CaviWipes™ Bleach	Compatible
Genius	CaviCide™	No information
Genius	CaviWipes™	No information
Genius	Envirocide™	No information
Genius	CaviCide1™	No information
Genius	CaviWipes1™	No information
Genius	CaviWipes™ 2.0	No information
Genius	CaviWipes™ HP	No information
Genius	CaviWipes™ Bleach	No information
Aixplorer SC6-1	CaviWipes1™	No information
Aixplorer SC6-1	CaviWipes™ 2.0	No information
Aixplorer SC6-1	CaviWipes™ HP	No information
Aixplorer SC6-1	CaviWipes™ Bleach	No information
Aixplorer SE 12-3	CaviCide™	No information
Aixplorer SE 12-3	CaviWipes™	No information
Aixplorer SE 12-3	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Aixplorer SE 12-3	CaviCide1™	No information
Aixplorer SE 12-3	CaviWipes1™	No information
Aixplorer SE 12-3	CaviWipes™ 2.0	No information
Aixplorer SE 12-3	CaviWipes™ HP	No information
Aixplorer SE 12-3	CaviWipes™ Bleach	No information
Aixplorer SEV 12-3	CaviCide™	No information
Aixplorer SEV 12-3	CaviWipes™	No information
Aixplorer SEV 12-3	Envirocide™	No information
Aixplorer SEV 12-3	CaviCide1™	No information
Aixplorer SEV 12-3	CaviWipes1™	No information
Aixplorer SEV 12-3	CaviWipes™ 2.0	No information
Aixplorer SEV 12-3	CaviWipes™ HP	No information
Aixplorer SEV 12-3	CaviWipes™ Bleach	No information
Aixplorer SEV 12-3 GENII	CaviCide™	No information
Aixplorer SEV 12-3 GENII	CaviWipes™	No information
Aixplorer SEV 12-3 GENII	Envirocide™	No information
Aixplorer SEV 12-3 GENII	CaviCide1™	No information
Aixplorer SEV 12-3 GENII	CaviWipes1™	No information
Aixplorer SEV 12-3 GENII	CaviWipes™ 2.0	No information
Aixplorer SEV 12-3 GENII	CaviWipes™ HP	No information
Aixplorer SEV 12-3 GENII	CaviWipes™ Bleach	No information
Aixplorer SL 10-2	CaviCide™	No information
Aixplorer SL 10-2	CaviWipes™	No information
Aixplorer SL 10-2	Envirocide™	No information
Aixplorer SL 10-2	CaviCide1™	No information
Aixplorer SL 10-2	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Aixplorer SL 10-2	CaviWipes™ 2.0	No information
Aixplorer SL 10-2	CaviWipes™ HP	No information
Aixplorer SL 10-2	CaviWipes™ Bleach	No information
Aixplorer SL 15-4	CaviCide™	No information
Aixplorer SL 15-4	CaviWipes™	No information
Aixplorer SL 15-4	Envirocide™	No information
Aixplorer SL 15-4	CaviCide1™	No information
Aixplorer SLV 16-5	CaviCide™	No information
Aixplorer SLV 16-5	CaviWipes™	No information
Aixplorer SLV 16-5	Envirocide™	No information
Aixplorer SLV 16-5	CaviCide1™	No information
Aixplorer SLV 16-5	CaviWipes1™	No information
Aixplorer SLV 16-5	CaviWipes™ 2.0	No information
Aixplorer SLV 16-5	CaviWipes™ HP	No information
Aixplorer SLV 16-5	CaviWipes™ Bleach	No information
Aixplorer SMC 12-3	CaviCide™	No information
Aixplorer SMC 12-3	CaviWipes™	No information
Aixplorer SMC 12-3	Envirocide™	No information
Aixplorer SMC 12-3	CaviCide1™	No information
Aixplorer SMC 12-3	CaviWipes1™	No information
Aixplorer SMC 12-3	CaviWipes™ 2.0	No information
Aixplorer SMC 12-3	CaviWipes™ HP	No information
Aixplorer SMC 12-3	CaviWipes™ Bleach	No information
Aixplorer XC6-1	CaviCide™	No information
Aixplorer XC6-1	CaviWipes™	No information
Aixplorer XC6-1	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Aixplorer XC6-1	CaviCide1™	No information
Aixplorer XC6-1	CaviWipes1™	No information
Aixplorer XC6-1	CaviWipes™ 2.0	No information
Aixplorer XC6-1	CaviWipes™ HP	No information
Aixplorer XC6-1	CaviWipes™ Bleach	No information
Aixplorer XP5-1	CaviCide™	No information
Aixplorer XP5-1	CaviWipes™	No information
Aixplorer XP5-1	Envirocide™	No information
Aixplorer XP5-1	CaviCide1™	No information
Aixplorer XP5-1	CaviWipes1™	No information
Aixplorer XP5-1	CaviWipes™ 2.0	No information
Aixplorer XP5-1	CaviWipes™ HP	No information
Aixplorer XP5-1	CaviWipes™ Bleach	No information
Aquilex	CaviCide™	No information
Aquilex	CaviWipes™	No information
Aquilex	Envirocide™	No information
Aquilex	CaviCide1™	No information
Aquilex	CaviWipes1™	No information
Aquilex	CaviWipes™ 2.0	No information
Aquilex	CaviWipes™ HP	No information
Aquilex	CaviWipes™ Bleach	No information
Imaging Equipment	CaviCide™	No information
Imaging Equipment	CaviWipes™	No information
Imaging Equipment	Envirocide™	No information
Imaging Equipment	CaviCide1™	No information
Imaging Equipment	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Imaging Equipment	CaviWipes™ 2.0	No information
Imaging Equipment	CaviWipes™ HP	Compatible
Imaging Equipment	CaviWipes™ Bleach	Compatible
Selenia	CaviCide™	No information
Selenia	CaviWipes™	No information
Selenia	Envirocide™	No information
Selenia	CaviCide1™	No information
Selenia	CaviWipes1™	No information
Selenia	CaviWipes™ 2.0	No information
Selenia	CaviWipes™ HP	Compatible
Selenia	CaviWipes™ Bleach	Compatible
Supersonic Imagine Aixplorer/exterior	CaviCide™	No information
Supersonic Imagine Aixplorer/exterior	CaviWipes™	No information
Supersonic Imagine Aixplorer/exterior	Envirocide™	No information
Supersonic Imagine Aixplorer/exterior	CaviCide1™	No information
Acessa ProVu	CaviCide™	No information
Acessa ProVu	CaviWipes™	No information
Acessa ProVu	Envirocide™	No information
Acessa ProVu	CaviCide1™	No information
Acessa ProVu	CaviWipes1™	No information
Acessa ProVu	CaviWipes™ 2.0	No information
Acessa ProVu	CaviWipes™ HP	No information
Acessa ProVu	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Affirm	CaviCide™	No information
Affirm	CaviWipes™	No information
Affirm	Envirocide™	No information
Affirm	CaviCide1™	No information
Affirm	CaviWipes1™	No information
Affirm	CaviWipes™ 2.0	No information
Affirm	CaviWipes™ HP	Compatible
Affirm	CaviWipes™ Bleach	Compatible
Aixplorer SC6-1	CaviCide™	No information
Aixplorer SC6-1	CaviWipes™	No information
Aixplorer SC6-1	Envirocide™	No information
Aixplorer SC6-1	CaviCide1™	No information
Aixplorer SL 15-4	CaviWipes1™	No information
Aixplorer SL 15-4	CaviWipes™ 2.0	No information
Aixplorer SL 15-4	CaviWipes™ HP	No information
Aixplorer SL 15-4	CaviWipes™ Bleach	No information
Aixplorer SL 18-5	CaviCide™	No information
Aixplorer SL 18-5	CaviWipes™	No information
Aixplorer SL 18-5	Envirocide™	No information
Aixplorer SL 18-5	CaviCide1™	No information
Aixplorer SL 18-5	CaviWipes1™	No information
Aixplorer SL 18-5	CaviWipes™ 2.0	No information
Aixplorer SL 18-5	CaviWipes™ HP	No information
Aixplorer SL 18-5	CaviWipes™ Bleach	No information
Aixplorer SLH 20-6	CaviCide™	No information
Aixplorer SLH 20-6	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Aixplorer SLH 20-6	Envirocide™	No information
Aixplorer SLH 20-6	CaviCide1™	No information
Aixplorer SLH 20-6	CaviWipes1™	No information
Aixplorer SLH 20-6	CaviWipes™ 2.0	No information
Aixplorer SLH 20-6	CaviWipes™ HP	No information
Aixplorer SLH 20-6	CaviWipes™ Bleach	No information
Supersonic Imagine Aixplorer/exterior	CaviWipes1™	No information
Supersonic Imagine Aixplorer/exterior	CaviWipes™ 2.0	No information
Supersonic Imagine Aixplorer/exterior	CaviWipes™ HP	No information
Supersonic Imagine Aixplorer/exterior	CaviWipes™ Bleach	No information
Supersonic Imagine Aixplorer/monitor	CaviCide™	No information
Supersonic Imagine Aixplorer/monitor	CaviWipes™	No information
Supersonic Imagine Aixplorer/monitor	Envirocide™	No information
Supersonic Imagine Aixplorer/monitor	CaviCide1™	No information
Supersonic Imagine Aixplorer/monitor	CaviWipes1™	No information
Supersonic Imagine Aixplorer/monitor	CaviWipes™ 2.0	No information
Supersonic Imagine Aixplorer/monitor	CaviWipes™ HP	No information
Supersonic Imagine Aixplorer/monitor	CaviWipes™ Bleach	No information

Equipment	Metrex Product	Compatibility Status
Supersonic Imagine Aixplorer/Scanning Housing	CaviCide™	Compatible
Supersonic Imagine Aixplorer/Scanning Housing	CaviWipes™	Compatible
Supersonic Imagine Aixplorer/Scanning Housing	Envirocide™	Compatible
Supersonic Imagine Aixplorer/Scanning Housing	CaviCide1™	No information
Supersonic Imagine Aixplorer/Scanning Housing	CaviWipes1™	No information
Supersonic Imagine Aixplorer/Scanning Housing	CaviWipes™ 2.0	No information
Supersonic Imagine Aixplorer/Scanning Housing	CaviWipes™ HP	Compatible
Supersonic Imagine Aixplorer/Scanning Housing	CaviWipes™ Bleach	Compatible
Thin Prep Review Scope	CaviCide™	No information
Thin Prep Review Scope	CaviWipes™	No information
Thin Prep Review Scope	Envirocide™	No information
Thin Prep Review Scope	CaviCide1™	No information
Thin Prep Review Scope	CaviWipes1™	No information
Thin Prep Review Scope	CaviWipes™ 2.0	No information
Thin Prep Review Scope	CaviWipes™ HP	No information
Thin Prep Review Scope	CaviWipes™ Bleach	No information
Trident	CaviCide™	No information
Trident	CaviWipes™	No information
Trident	Envirocide™	No information
Trident	CaviCide1™	No information
Trident	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Trident	CaviWipes™ 2.0	No information
Trident	CaviWipes™ HP	Compatible
Trident	CaviWipes™ Bleach	Compatible
Novasure Model 10RF	CaviCide™	No information
Novasure Model 10RF	CaviWipes™	No information
Novasure Model 10RF	Envirocide™	No information
Novasure Model 10RF	CaviCide1™	No information
Novasure Model 10RF	CaviWipes1™	No information
Novasure Model 10RF	CaviWipes™ 2.0	No information
Novasure Model 10RF	CaviWipes™ HP	No information
Novasure Model 10RF	CaviWipes™ Bleach	No information

Masimo		
Equipment	Metrex Product	Compatibility Status
ECG Sensors	CaviCide™	No information
ECG Sensors	CaviWipes™	No information
ECG Sensors	Envirocide™	No information
ECG Sensors	CaviCide1™	No information
ECG Sensors	CaviWipes1™	No information
ECG Sensors	CaviWipes™ 2.0	No information
ECG Sensors	CaviWipes™ HP	No information
ECG Sensors	CaviWipes™ Bleach	No information
Radical-7	CaviCide™	Compatible
Radical-7	CaviWipes™	Compatible
Radical-7	Envirocide™	Compatible
Radical-7	CaviCide1™	Compatible
Radical-7	CaviWipes1™	Compatible

Equipment	Metrex Product	Compatibility Status
Radical-7	CaviWipes™ 2.0	Compatible
Radical-7	CaviWipes™ HP	No information
Radical-7	CaviWipes™ Bleach	Compatible