

Material Compatibility Status

DentaleEZ		
Equipment	Metrex Product	Compatibility Status
Plastic & coated Metal surfaces of dental chairs	CaviCide™	Compatible
Plastic & coated Metal surfaces of dental chairs	CaviWipes™	Compatible
Plastic & coated Metal surfaces of dental chairs	Envirocide™	Compatible
Plastic & coated Metal surfaces of dental chairs	CaviCide1™	No information
Plastic & coated Metal surfaces of dental chairs	CaviWipes1™	No information
Plastic & coated Metal surfaces of dental chairs	CaviWipes™ 2.0	No information
Plastic & coated Metal surfaces of dental chairs	CaviWipes™ HP	No information
Plastic & coated Metal surfaces of dental chairs	CaviWipes™ Bleach	No information
NuSmplicity Chair	CaviCide™	Compatible
NuSmplicity Chair	CaviWipes™	Compatible
NuSmplicity Chair	Envirocide™	Compatible
NuSmplicity Chair	CaviCide1™	No information
NuSmplicity Chair	CaviWipes1™	No information
NuSmplicity Chair	CaviWipes™ 2.0	No information
NuSmplicity Chair	CaviWipes™ HP	No information
NuSmplicity Chair	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
Upholstery	CaviCide™	Compatible
Upholstery	CaviWipes™	Compatible
Upholstery	Envirocide™	Compatible
Upholstery	CaviCide1™	No information
Upholstery	CaviWipes1™	No information
Upholstery	CaviWipes™ 2.0	No information
Upholstery	CaviWipes™ HP	No information
Upholstery	CaviWipes™ Bleach	Compatible

Furniss Corporation		
Equipment	Metrex Product	Compatibility Status
Phoenix Model 1850 Knee CPM	CaviCide™	No information
Phoenix Model 1850 Knee CPM	CaviWipes™	No information
Phoenix Model 1850 Knee CPM	Envirocide™	No information
Phoenix Model 1850 Knee CPM	CaviCide1™	No information
Phoenix Model 1850 Knee CPM	CaviWipes1™	No information
Phoenix Model 1850 Knee CPM	CaviWipes™ 2.0	No information
Phoenix Model 1850 Knee CPM	CaviWipes™ HP	No information
Phoenix Model 1850 Knee CPM	CaviWipes™ Bleach	No information

Implant Direct		
-----------------------	--	--

Equipment	Metrex Product	Compatibility Status
Drill Stop Kit	CaviCide™	No information
Drill Stop Kit	CaviWipes™	No information
Drill Stop Kit	Envirocide™	No information
Drill Stop Kit	CaviCide1™	No information
Drill Stop Kit	CaviWipes1™	No information
Drill Stop Kit	CaviWipes™ 2.0	No information
Drill Stop Kit	CaviWipes™ HP	No information
Drill Stop Kit	CaviWipes™ Bleach	No information
Mini Sterilizable Tray	CaviCide™	No information
Mini Sterilizable Tray	CaviWipes™	No information
Mini Sterilizable Tray	Envirocide™	No information
Mini Sterilizable Tray	CaviCide1™	No information
Mini Sterilizable Tray	CaviWipes1™	No information
Mini Sterilizable Tray	CaviWipes™ 2.0	No information
Mini Sterilizable Tray	CaviWipes™ HP	No information
Mini Sterilizable Tray	CaviWipes™ Bleach	No information
Standard Sterizable Trays	CaviCide™	No information
Standard Sterizable Trays	CaviWipes™	No information
Standard Sterizable Trays	Envirocide™	No information
Standard Sterizable Trays	CaviCide1™	No information
Standard Sterizable Trays	CaviWipes1™	No information
Standard Sterizable Trays	CaviWipes™ 2.0	No information
Standard Sterizable Trays	CaviWipes™ HP	No information
Standard Sterizable Trays	CaviWipes™ Bleach	No information