

Material Compatibility Status

Nova Biomedical		
Equipment	Metrex Product	Compatibility Status
StatStrip Express	CaviCide™	No information
StatStrip Express	CaviWipes™	No information
StatStrip Express	Envirocide™	No information
StatStrip Express	CaviCide1™	No information
StatStrip Express	CaviWipes1™	No information
StatStrip Express	CaviWipes™ 2.0	No information
StatStrip Express	CaviWipes™ HP	No information
StatStrip Express	CaviWipes™ Bleach	Compatible