

Material Compatibility Status

GE Healthcare		
Equipment	Metrex Product	Compatibility Status
3.0T Signa HDxt	CaviCide™	Not recommended for use
3.0T Signa HDxt	CaviWipes™	Not recommended for use
3.0T Signa HDxt	Envirocide™	Not recommended for use
3.0T Signa HDxt	CaviCide1™	Not recommended for use
3.0T Signa HDxt	CaviWipes1™	Not recommended for use
3.0T Signa HDxt	CaviWipes™ 2.0	Not recommended for use
3.0T Signa HDxt	CaviWipes™ HP	No information
3.0T Signa HDxt	CaviWipes™ Bleach	Compatible
3Sc-RS	CaviCide™	Compatible
3Sc-RS	CaviWipes™	Compatible
3Sc-RS	Envirocide™	Compatible
3Sc-RS	CaviCide1™	Compatible
3Sc-RS	CaviWipes1™	Compatible
3Sc-RS	CaviWipes™ 2.0	No information
3Sc-RS	CaviWipes™ HP	No information
3Sc-RS	CaviWipes™ Bleach	No information