

Material Compatibility Status

GE Healthcare		
Equipment	Metrex Product	Compatibility Status
Discovery MR750 3.0T	CaviCide™	Not recommended for use
Discovery MR750 3.0T	CaviWipes™	Not recommended for use
Discovery MR750 3.0T	Envirocide™	Not recommended for use
Discovery MR750 3.0T	CaviCide1™	Not recommended for use
Discovery MR750 3.0T	CaviWipes1™	Not recommended for use
Discovery MR750 3.0T	CaviWipes™ 2.0	Not recommended for use
Discovery MR750 3.0T	CaviWipes™ HP	No information
Discovery MR750 3.0T	CaviWipes™ Bleach	Compatible
RIC6=12-D	CaviCide™	Compatible
RIC6=12-D	CaviWipes™	Compatible
RIC6=12-D	Envirocide™	Compatible
RIC6=12-D	CaviCide1™	Compatible
RIC6=12-D	CaviWipes1™	Compatible
RIC6=12-D	CaviWipes™ 2.0	No information
RIC6=12-D	CaviWipes™ HP	No information
RIC6=12-D	CaviWipes™ Bleach	No information