

# Material Compatibility Status

| GE Healthcare |                   |                         |
|---------------|-------------------|-------------------------|
| Equipment     | Metrex Product    | Compatibility Status    |
| Optima MR360  | CaviWipes1™       | Not recommended for use |
| Optima MR360  | CaviWipes™ 2.0    | Not recommended for use |
| Optima MR360  | CaviWipes™ HP     | No information          |
| Optima MR360  | CaviWipes™ Bleach | Compatible              |
| Optima MR360  | CaviCide™         | Not recommended for use |
| Optima MR360  | CaviWipes™        | Not recommended for use |
| Optima MR360  | Envirocide™       | Not recommended for use |
| Optima MR360  | CaviCide1™        | Not recommended for use |