

Material Compatibility Status

GE Healthcare		
Equipment	Metrex Product	Compatibility Status
Dash 4000	CaviCide™	Not recommended for use
Dash 4000	CaviWipes™	Not recommended for use
Dash 4000	Envirocide™	Not recommended for use
Dash 4000	CaviCide1™	Not recommended for use
Dash 4000	CaviWipes1™	Not recommended for use
Dash 4000	CaviWipes™ 2.0	Not recommended for use
Dash 4000	CaviWipes™ HP	No information
Dash 4000	CaviWipes™ Bleach	Compatible