

Material Compatibility Status

GE Healthcare		
Equipment	Metrex Product	Compatibility Status
C2-6B-D	CaviCide™	Compatible
C2-6B-D	CaviWipes™	Compatible
C2-6B-D	Envirocide™	Compatible
C2-6B-D	CaviCide1™	No information
C2-6B-D	CaviWipes1™	No information
C2-6B-D	CaviWipes™ 2.0	No information
C2-6B-D	CaviWipes™ HP	No information
C2-6B-D	CaviWipes™ Bleach	No information
Vivid T8	CaviCide™	Compatible
Vivid T8	CaviWipes™	Compatible
Vivid T8	Envirocide™	Compatible
Vivid T8	CaviCide1™	No information
Vivid T8	CaviWipes1™	No information
Vivid T8	CaviWipes™ 2.0	No information
Vivid T8	CaviWipes™ HP	No information
Vivid T8	CaviWipes™ Bleach	No information